NGeneAnalySys™

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| Category | Database | Version | Release | Description | http |
|-------------------------|------------------------------------|---------|----------|--|------|
| Disease Database | ClinVar | 201601 | Jan.2016 | Database of assertions about the clinical significance and phenotype relationship of human variations. | |
| Disease Database (BRCA) | ARUP-BRCA Mutation Database (ARUP) | . | Sep.2017 | Through collaborations with the University of Utah Huntsman Cancer Institute (HCI), and with the WHO International Agency for Research on Cancer (IARC), the University of Utah Department of Pathology and ARUP Laboratories is pleased to host the BRCA mutation database. | |
| Disease Database (BRCA) | BRCA Exchange | v3.33.0 | Jan.2017 | curated data platform that provides information on variants of specific genes: BRCA1 and BRCA2. | |
| Disease Database (BRCA) | BRCA Share | . | Jun.2017 | BRCA Share (formerly UMD-Universal Mutation Database) | |

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| Software | License | Version | Release | Description | http |
|----------------|----------------------------|------------|----------|--|------|
| BWA | GPLv3, MIT | v0.7.10 | Jul.2014 | mapping low-divergent sequences against a large reference genome | |
| FLASH | GPLv3 | v1.2.11 | Aug.2014 | merge paired-end reads | |
| FastQC | GPLv3 | v0.11.3 | May.2015 | quality control tool for high throughput sequence data | |
| FastQValidator | the University of Michigan | v0.1.1 | Oct.2012 | validation the format of FASTQ file | |
| GATK | MIT | v2.3 Lite | Dec.2012 | wide variety of tools with a primary focus on variant discovery and genotyping | |
| GMAP | Thomas D. Wu | 2016-11-07 | Nov.2016 | A Genomic Mapping and Alignment Program for mRNA and EST Sequences | |

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The screenshot shows the 'New Analysis Request' page from NGeneAnalySys. The interface includes a header with the application name and a navigation link. Below the header, there are several input fields and dropdown menus. Callouts numbered 1 through 5 highlight specific elements:

- 1**: A red box surrounds the 'Run' field, which contains the placeholder 'Enter the Run Name'.
- 2**: A red box surrounds the 'Sequencer' dropdown menu, which lists 'Illumina MiSeqDX', 'Illumina MiSeq', and 'Illumina NextSeq 550Dx'.
- 3**: A red box surrounds the 'FASTQ 파일 선택' (FASTQ File Selection) section, which includes three radio buttons: 'Local Fastq Files' (selected), 'Server Fastq Files', and 'Server Run Folder'.
- 4**: A large red box surrounds the '검체 정보' (Sample Information) table, which contains four rows: Sample (HD795-2, HD795-5), Panel (BRCAaccuTest PLUS v1, BRCAaccuTest PLUS v1), Source (BLOOD, BLOOD), and Disease (HBOC, HBOC).
- 5**: A red box surrounds the 'SUBMIT' button at the bottom right of the form.

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
Q Queued R Running C Complete F Failed

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The screenshot shows the NGeneAnalySys software interface. The top navigation bar has tabs: DASHBOARD, RESULTS (which is highlighted with a yellow circle and red box), SAMPLES, and MANAGER. On the far right, there's an administrator dropdown and a user icon. The main content area is titled "RESULTS". It displays three tables of analysis results, each with a header row and multiple rows of data. The first table is for run 2018-04-30:09:24:50, the second for 2018-04-30:09:24:38, and the third for 2018-04-30:09:23:46. Each table has columns for Run, Sequencer, Start Date, Finished Date, Sample, Panel, Status, Variant, and QC. The search panel on the right allows filtering by DATE (two date pickers), PANEL (a dropdown with options like BRCAaccuTest v2 and BRCAaccuTest PLUS v1), and SAMPLE (a text input with a date range from 2018-04-01 to 2018-04-30). A version number v 1.3 (Build Date 20180412_3) is at the bottom right.

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The screenshot shows the NGeneAnalySys software interface. The main window displays search results for sequencing runs. The search panel on the right is highlighted with a red box and numbered 1 through 4. Number 1 points to the search icon in the header. Number 2 points to the date input fields. Number 3 points to the search button. Number 4 points to the search results table.

| Run | Sequencer | Start Date | Finished Date |
|---------------------|---------------------|-----------------------------------|---------------|
| 2018-04-30:09:24:50 | Illumina MiSeqDx | 2018/04/30 | 2018/04/30 |
| Mix-12-5-each | BRCAccuTest PLUS v1 | C: 0 P: 0 LP: 0 US: 0 LB: 0 B: 12 | PASS |
| Mix-2-each | BRCAccuTest PLUS v1 | C: 0 P: 0 LP: 0 US: 0 LB: 0 B: 12 | PASS |
| Mix-5-each | BRCAccuTest PLUS v1 | C: 0 P: 0 LP: 0 US: 0 LB: 0 B: 12 | PASS |
| 2018-04-30:09:24:38 | Illumina MiSeqDx | 2018/04/30 | 2018/04/30 |
| HD795-2 | BRCAccuTest PLUS v1 | C: 0 P: 0 LP: 0 US: 0 LB: 0 B: 12 | PASS |
| HD795-5 | BRCAccuTest PLUS v1 | C: 0 P: 2 LP: 0 US: 1 LB: 1 B: 9 | PASS |
| 2018-04-30:09:23:46 | Illumina MiSeqDx | 2018/04/30 | 2018/04/30 |

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The screenshot shows the NGeneAnalySys software interface. At the top, there are tabs: DASHBOARD, RESULTS (highlighted in green), SAMPLES, and MANAGER. On the right, there's a sidebar for 'Administrator' and 'NGeneBio Innovative & Next Generation Biotechnology'. The main area displays a results table with the following data:

| Run | Sequencer | Start Date | Finished Date | |
|---------------------|----------------------|----------------|--|------|
| 2018-04-30:09:24:50 | Illumina MiSeqDx | 2018/04/30 | 2018/04/30 | |
| Sample | Panel | Status | Variant | QC |
| Mix-12-5-each | BRCAaccuTest PLUS v1 | C | P 0 LP 0 US 0 LB 0 B 12 | PASS |
| Mix-2-each | BRCAaccuTest PLUS v1 | C | P 0 LP 0 US 0 LB 0 B 12 | PASS |
| Mix-5-each | BRCAaccuTest PLUS v1 | C | P 0 LP 0 US 0 LB 0 B 12 | PASS |

Below this table are two smaller tables:

| Sample | Panel | Status | Variant | QC |
|---------|----------------------|----------------|--|------|
| HD795-2 | BRCAaccuTest PLUS v1 | C | P 0 LP 0 US 0 LB 0 B 12 | PASS |
| HD795-5 | BRCAaccuTest PLUS v1 | C | P 2 LP 0 US 1 LB 1 B 9 | PASS |

At the bottom, there are search filters for DATE, PANEL, and DATE range, along with a version note: v 1.3 (Build Date 20180412_3).

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
The screenshot shows the NGeneAnalySys™ software interface with the following details:

- Header:** NGeneAnalySys™, DASHBOARD, RESULTS, SAMPLES (highlighted), MANAGER, Administrotor.
- Top Bar:** NA14805-1, RUN : BRCA+ Test, PANEL : BRCAaccuTest PLUS v1, DISEASE : HBOC, SEQUENCER : Illumina MiSeqDx.
- Section 1:** OVERVIEW (highlighted with a red box and circled with a yellow circle).
- Section 2:** DETECTED VARIANT SUMMARY (highlighted with a red box and circled with a yellow circle).
 - P: Variants : 1, Gene : 1
 - LP: Variants : 0, Gene : 0
 - US: Variants : 0, Gene : 0
 - LB: Variants : 1, Gene : 1
 - B: Variants : 13, Gene : 2
- Section 3:** Pathogenicity, Gene, Position, Variant, Transcript, NT Change, AA Change (highlighted with a red box and circled with a yellow circle).

| Pathogenicity | Gene | Position | Variant | Transcript | NT Change | AA Change |
|---------------|-------|----------|----------------------|-------------|----------------------|---------------------------|
| P | BRCA2 | 32900700 | NM_000059.3:c.581G>A | NM_000059.3 | NM_000059.3:c.581G>A | NP_000050.2:p.(Trp194Ter) |
- Section 4:** QC METRICS (highlighted with a red box and circled with a yellow circle).

| ROI Coverage | Coverage Uniformity |
|--------------|---------------------|
| PASS | PASS |

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| 분류 | 아이콘 | 설명 |
|------------------------|-----|--|
| Pathogenic | P | Causative for clinical phenotype. |
| Likely Pathogenic | LP | Probably causative for clinical phenotype. |
| Uncertain Significance | US | Unknown effect on clinical phenotype. |
| Likely Benign | LB | Probably not causative for a clinical phenotype. |
| Benign | B | Not causative for a clinical phenotype. |

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| SNV/Indel : BRCA1 synonymous_variant NM_007294.3 c.2082C>T p.(S694=) | | | | | | | | | |
|--|---------------|---------|--------|-------|-------------|------|--------------------|-----------------------|-----------------------|
| Prediction | Pathogenicity | Warning | Report | Gene | Transcript | Type | Consequence | NT Change | NT Change (BIC) |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.2082C>T | NM_007294.3:c.2201... |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.2311T>C | NM_007294.3:c.2430... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.2612C>T | NM_007294.3:c.2731... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.3113A>G | NM_007294.3:c.3232... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.3548A>G | NM_007294.3:c.3667... |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.4308T>C | NM_007294.3:c.4427... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.4837A>G | NM_007294.3:c.4956... |
| LB | B | | R | BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.5744C>T | NM_000059.3:c.5972... |
| P | | | | BRCA2 | NM_000059.3 | SNP | stop_gained | NM_000059.3:c.581G>A | NM_000059.3:c.809... |
| R | | | | BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.720T>C | NM_000059.3:c.7625... |
| | | | | BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.7470... | NM_000059.3:c.7470... |

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| SNV/Indel : BRCA2 stop_gained NM_000059.3 c.581G>A p.(W194Ter) | | | | | | | | | |
|--|---------------|---------|------------------|-----------------------|-----------------------|------|--------------------|-----------------------|-----------------------|
| Prediction | Pathogenicity | Warning | Report | Gene | Transcript | Type | Consequence | NT Change (BIC) | AA Change |
| P | | | | BRCA2 | NM_000059.3 | SNP | stop_gained | NM_000059.3:c.809... | NM_000059.3:c.809... |
| LB | B | | R | BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.5972... | NM_000059.3:c.5972... |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.2430... | NM_007294.3:c.2430... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.2612C>T | NM_007294.3:c.2731... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.3113A>G | NM_007294.3:c.3232... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.3548A>G | NM_007294.3:c.3667... |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.4308T>C | NM_007294.3:c.4427... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.4837A>G | NM_007294.3:c.4956... |
| BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.5744C>T | NM_000059.3:c.5972... | | | | |
| BRCA2 | NM_000059.3 | SNP | stop_gained | NM_000059.3:c.581G>A | NM_000059.3:c.809... | | | | |
| BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.720T>C | NM_000059.3:c.7625... | | | | |
| BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.7470... | NM_000059.3:c.7470... | | | | |

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| SNV/Indel : BRCA2 stop_gained NM_000059.3 c.581G>A p.(W194Ter) | | | | | | | | | | |
|--|---------------|---------|--------|-------|-------------|------|--------------------|-----------------------|-----------------------|------------|
| Prediction | Pathogenicity | Warning | Report | Gene | Transcript | Type | Consequence | NT Change | NT Change (BIC) | AA Change |
| P | | | | BRCA2 | NM_000059.3 | SNP | stop_gained | NM_000059.3:c.581G>A | NM_000059.3:c.809... | NP_0000... |
| LB | B | R | | BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.574C>T | NM_000059.3:c.597... | NP_0000... |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.2311T>C | NM_007294.3:c.2430... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.2612C>T | NM_007294.3:c.2731... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.3113A>G | NM_007294.3:c.3232... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.3548A>G | NM_007294.3:c.3667... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.4308T>C | NM_007294.3:c.4427... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.4837A>G | NM_007294.3:c.4956... | NP_0092... |
| B | | | | BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.7397C= | NM_000059.3:c.7625... | NP_0000... |
| | | | | BRCA2 | NM_000059.3 | SNP | synonymous_variant | NM_000059.3:c.7747A>C | NM_000059.3:c.7747A>C | NP_0000... |

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
NGeneAnalySys > Variant Filter

Filter 1

Variant Consequence Population Frequency

Prediction Tier1 Tier2 Tier3 Tier4

Pathogenicity P(Pathogenic) LP(Likely Pathogenic) US(Uncertain Significance) LB(Likely Benign) B(Benign)

ClinVar P(Pathogenic) LP(Likely Pathogenic) US(Uncertain Significance) LB(Likely Benign) B(Benign)

Variant Type SNP Ins(Insertion) Del(Deletion)
 MNP(Multiple-nucleotide polymorphism) Complex(Multiple-nucleotides and an InDel)

Chromosome chr7, chr9

Gene FLT3, CEBPA

dbSNP Has dbSNP ID

COSMIC Has COSMIC ID


COSMIC Occurrence Tissue selection

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


The screenshot shows the 'Variant Filter' dialog box from the NGeneAnalySys software. The 'Filter' dropdown menu is set to 'TEST' and has a red border around it. At the bottom right of the dialog, there are two buttons: 'Cancel' and 'Save Filter', with 'Save Filter' also having a red border around it. The dialog contains several sections for filtering variants:


- Variant**: A tab selected at the top.
- Consequence**: A tab next to Variant.
- Population Frequency**: A tab next to Consequence.
- Prediction**: Options for Tier1, Tier2, Tier3, and Tier4.
- Pathogenicity**: Options for P(Pathogenic), LP(Likely Pathogenic), US(Uncertain Significance), LB(Likely Benign), and B(Benign).
- ClinVar**: Options for P(Pathogenic), LP(Likely Pathogenic), US(Uncertain Significance), LB(Likely Benign), and B(Benign).
- Variant Type**: Options for SNP, Ins(Insertion), Del(Deletion), MNP(Multiple-nucleotide polymorphism), and Complex(Multiple-nucleotides and an InDel).
- Chromosome**: An input field containing 'chr7, chr9'.
- Gene**: An input field containing 'FLT3, CEBPA'.
- dbSNP**: A checkbox labeled 'Has dbSNP ID'.
- COSMIC**: A checkbox labeled 'Has COSMIC ID'.
- COSMIC Occurrence**: A dropdown menu set to 'Tissue selection'.
- Allele Fraction**: Input fields for '30' and '50' followed by a '%' symbol.

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

The screenshot shows a user interface for analyzing SNV/Indel variants. At the top left, there's a dropdown menu labeled "Clinical Significant". Below it, two main sections are displayed within a red-bordered box:

- Predicted Pathogenicity:** Shows four sets of colored boxes representing predictions from different tools: P (grey), LP (light blue), US (white), LB (dark blue), and B (dark grey). The first set is all grey (P, LP, US, LB, B). The second set has LP and US white, while P, LB, and B are dark. The third set has LP and US white, while P, LB, and B are dark. The fourth set has LP and US white, while P, LB, and B are dark.
- Change Pathogenicity:** Shows four sets of colored boxes representing changes in pathogenicity from different tools: P (grey), LP (light blue), US (white), LB (dark blue), and B (dark grey). The first set is all grey (P, LP, US, LB, B). The second set has LP and US light blue, while P, LB, and B are dark. The third set has LP and US white, while P, LB, and B are dark. The fourth set has LP and US white, while P, LB, and B are dark.


At the top right, there's a yellow circle with the number "3" and a button labeled "Add to Report" with a checked checkbox.

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


| SNV/Indel : BRCA2 stop_gained NM_000059.3 c.581G>A p.(W194Ter) | | | | | | | | | | |
|--|---------------|---------|--------|-------|-------------|------|--------------------|-----------------------|-----------------------|------------|
| Prediction | Pathogenicity | Warning | Report | Gene | Transcript | Type | Consequence | NT Change | NT Change (BIC) | AA Change |
| P | | | | BRCA2 | NM_000059.3 | SNP | stop_gained | NM_000059.3:c.581G>A | NM_000059.3:c.809... | NP_0000... |
| LB | | | | BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.5744C>T | NM_000059.3:c.5972... | NP_0000... |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.2311T>C | NM_007294.3:c.2430... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.2612C>T | NM_007294.3:c.2731... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.3113A>G | NM_007294.3:c.3232... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.3548A>G | NM_007294.3:c.3667... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | synonymous_variant | NM_007294.3:c.4308T>C | NM_007294.3:c.4427... | NP_0092... |
| B | | | | BRCA1 | NM_007294.3 | SNP | missense_variant | NM_007294.3:c.4837A>G | NM_007294.3:c.4956... | NP_0092... |
| B | | | | BRCA2 | NM_000059.3 | SNP | missense_variant | NM_000059.3:c.7397C= | NM_000059.3:c.7625... | NP_0000... |
| B | | | | BRCA2 | NM_000059.3 | SNP | synonymous_variant | NM_000059.3:c.7347A>C | NM_000059.3:c.7470... | NP_0000... |

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



The screenshot shows an interpretation log entry with the following details:

| Date | Type | User | Previous | Current | Comment |
|---------------------|------|-------|----------|---------|----------------------|
| 2018-05-16 11:32:43 | Tier | admin | T1 | T1 | therapeutic, A, Test |

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| Clinical Report | | | | | |
|--|----------------------------|--|--------------|------|----------------|
| Test Information | Client Information | | | | |
| Test No. | Ordering Organization | | | | |
| Test Name | Ordering Physician | | | | |
| Specimen | Patient | | | | |
| Specimen Type | BLOOD | | | | |
| Draw Date | Date of Birth | | | | |
| Received Date | Patient ID | | | | |
| Gender | | | | | |
| Summary | | | | | |
| Positive clinically significant mutation identified. | | | | | |
| BRCA1, NM_000293.3:c.514C>T(NP_000292.1:p.Thr171Met), heterozygote | | | | | |
| Result/Conclusion | | | | | |
| Details, BRCA1 | | | | | |
| NT Change (BRCA Format) | NT Change (HGVS Format) | AA Change* | Zygosity | Type | Classification |
| NM_007294.3:c.2 | NM_007294.3:c.2 | NP_000225.1:p. Gln117C>T (Leu117Ter) | heterozygote | SNP | B |
| ATG>TC | ATG>TC | NP_000225.1:p. Leu117Ter | heterozygote | SNP | B |
| TTC>CT | TTC>CT | NP_000225.1:p. (Ter117) | heterozygote | SNP | B |
| <small>Legend: B=benign, C=suspicious, H=harmful, N=negligible, U=uncertain, X=indeterminate, S=suspected pathogenic variant, P=pathogenic variant, V=variant of unknown significance, I=intermediate, A=likely benign, L=likely pathogenic, D=disease causing, M=mild, M=moderate, M=severe, M=very severe, M=extreme. *AA change is based on the reference sequence NM_000293.3. The AA change is bolded if it is a frameshift mutation.</small> | | | | | |
| Test Information | | | | | |
| 1. Specimen | BLOOD | | | | |
| 2. Clinical information | | | | | |
| 3. Test method | | | | | |
| Details, BRCA2 | | | | | |
| NT Change (BRCA Format) | NT Change (HGVS Format) | AA Change* | Zygosity | Type | Classification |
| NM_000293.3:c.8 | NM_000293.3:c.8 | NP_000220.2:p. Ile104>A | heterozygote | SNP | P |
| ATG>AG | ATG>AG | NP_000220.2:p. Ile104Ter | heterozygote | SNP | EB |
| NM_000293.3:c.5 | NM_000293.3:c.5 | NP_000220.2:p. Ile104Ter | heterozygote | SNP | EB |
| ATG>AG | ATG>AG | NP_000220.2:p. Ile104Ter | heterozygote | SNP | EB |
| NM_000293.3:c.13 | NM_000293.3:c.13 | NP_000220.2:p. Asp446>G | heterozygote | SNP | B |
| ATG>CG | ATG>CG | NP_000220.2:p. Asp446Ter | heterozygote | SNP | B |
| NM_000293.3:c.6 | NM_000293.3:c.6 | NP_000220.2:p. Ile142>V | heterozygote | SNP | B |
| ATG>GT | ATG>GT | NP_000220.2:p. Ile142Ter | heterozygote | SNP | B |
| NM_000293.3:c.4 | NM_000293.3:c.4 | NP_000220.2:p. Leu152I | heterozygote | SNP | B |
| ATG>CG | ATG>CG | NP_000220.2:p. Leu152Ter | heterozygote | SNP | B |
| NM_000293.3:c.3 | NM_000293.3:c.3 | NP_000220.2:p. Ile113D | heterozygote | SNP | B |
| ATG>CG | ATG>CG | NP_000220.2:p. (Ter113) | heterozygote | SNP | B |
| <small>Legend: B=benign, C=suspicious, H=harmful, N=negligible, U=uncertain, X=indeterminate, S=suspected pathogenic variant, P=pathogenic variant, V=variant of unknown significance, I=intermediate, A=likely benign, L=likely pathogenic, D=disease causing, M=mild, M=moderate, M=severe, M=very severe, M=extreme. *AA change is based on the reference sequence NM_000293.3. The AA change is bolded if it is a frameshift mutation.</small> | | | | | |
| Test Information | | | | | |
| 1. Specimen | BLOOD | | | | |
| 2. Clinical information | | | | | |
| 3. Test method | | | | | |

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The screenshot shows the NGeneAnalySys™ software interface. The top navigation bar includes DASHBOARD, RESULTS, SAMPLES, and MANAGER. The current view is the OVERVIEW tab. The main content area displays a summary of detected variants across four tiers:

- Tier I:** Variants: 2, Gene: 1
- Tier II:** Variants: 5, Gene: 4
- Tier III:** Variants: 52, Gene: 38
- Tier IV:** Variants: 154, Gene: 66

Below this summary is a detailed table of variants:

| Tier | Gene | Position | Variant | Transcript | NT Change | AA Change | VAF(Depth) | Evidence |
|---------|-------|----------|----------------|-------------|----------------|-------------|------------|----------|
| Tier I | RUNX1 | 36206711 | c.801G>A | NM_001754 | c.801G>A | p.Met267Ile | 53.99 | A, C |
| Tier I | RUNX1 | 36259345 | c.146C>T | NM_001754 | c.146C>T | p.Pro49Leu | 55.22 | A, C |
| Tier II | ETV6 | 12006451 | c.427_428delCA | NM_001987 | c.427_428delCA | p.Gln143fs | 6.50 | C |
| Tier II | BCOR | 20923086 | c.3621dninA | NM_00123385 | c.3621dninA | n.Gln1208fs | 82.32 | C |

At the bottom, the QC Metrics section shows the following results:

| Duplicated Reads | Mapped Base | On Target | On Target Coverage | Q30 Trimmed Base | ROI Coverage | Total Base |
|------------------|-------------|-----------|--------------------|------------------|--------------|------------|
| PASS | PASS | PASS | PASS | PASS | PASS | PASS |

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The screenshot shows a table of SNV/Indel variants for gene ABL1. The columns include Prediction, Tier, Warning, Low Conf., Report, Gene, Transcript, Type, Consequence, NT Change, and AA Change (Single). The 'Gene' column header is highlighted with a red box. Red warning icons are visible in the 'Tier' and 'Warning' columns for several rows.

| SNV/Indel : ABL1 missense_variant NM_007313 c.827A>G p.Y276C | | | | | | | | | | |
|--|------|---------|-------------|--------|------|--------------|------|--------------------|-----------|--------------------|
| Prediction | Tier | Warning | Low Conf... | Report | Gene | Transcript | Type | Consequence | NT Change | AA Change (Single) |
| T3 | R | | | | ABL1 | NM_007313 | SNP | missense_variant | c.827A>G | p.Y276C |
| T3 | R | | | | ABL2 | NM_001136000 | SNP | missense_variant | c.35G>C | p.S12T |
| T4 | | | | | ABL2 | NM_007314 | SNP | synonymous_variant | c.1482C>T | p.V494V |
| T4 | | ⚠️ | | | ABL2 | NM_007314 | SNP | synonymous_variant | c.2670T>G | p.G890G |
| T4 | | ⚠️ | | | AKT1 | NM_001014431 | SNP | synonymous_variant | c.726G>A | p.E242E |
| T4 | | | | | ALK | NM_004304 | SNP | missense_variant | c.4472A>G | p.K1491R |
| T4 | | | | | ALK | NM_004304 | SNP | missense_variant | c.4381A>G | p.I1461V |
| T4 | | | | | ALK | NM_004304 | SNP | synonymous_variant | c.2535T>C | p.G845G |
| T4 | | | | | ALK | NM_004304 | SNP | missense_variant | c.4587C>G | p.D1529E |
| T4 | | | | | ALK | NM_004304 | SNP | synonymous_variant | c.777E>C | p.C205C |

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The screenshot shows a table of SNV/Indel variants for gene RUNX1. The 'AA Change' column is expanded, indicated by a red box and a red arrow pointing to the plus sign icon in the header. The expanded list includes AA Change (Single), COSMIC ID, dbSNP ID (highlighted in blue), dbSNP Common ID, Exon, Chr, Start Position, Ref, Alt, Fraction, and Depth.

| SNV/Indel : RUNX1 missense_variant NM_001754 c.801G>A p.M267I | | | | | | | | | |
|---|------|---------|-------------|--------|-------|--------------|------|----------------------|--------------------|
| Prediction | Tier | Warning | Low Conf... | Report | Gene | Transcript | Type | AA Change | AA Change (Single) |
| T1 | R | | | | RUNX1 | NM_001754 | SNP | ✓ AA Change (Single) | p.M267I |
| T1 | R | | | | RUNX1 | NM_001754 | SNP | ✓ COSMIC ID | p.P49L |
| T2 | T1 | | | | KRAS | NM_033360 | SNP | ✓ dbSNP ID | p.G13D |
| T2 | | | | | BCOR | NM_00123385 | Ins | dbSNP Common ID | p.Q1208fs |
| T2 | | | | | NF1 | NM_001042492 | Del | ✓ Exon | p.P678fs |
| T2 | | | | | NF1 | NM_001042492 | Ins | ✓ Chr | p.I679fs |
| T2 | | | | | ETV6 | NM_001987 | Del | Start Position | p.Q143fs |
| T2 | | | | | KMT2D | NM_003482 | Ins | Ref | p.R2443fs |
| T3 | | | | | IL2RB | NM_000878 | Del | Alt | p.P489fs |
| T3 | | | | | NF1 | NM_001042492 | CNV | Fraction | p.D280T |
| T3 | | | | | | | | Depth | |

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The screenshot shows the NGeneAnalySys software interface with the 'SAMPLES' tab selected. On the left, there is a 'Filter' button. Three numbered circles (1, 2, 3) point to specific areas: 1 points to the 'Select Filter' dropdown menu, which is set to 'TEST'; 2 points to the 'Edit Filter' button; and 3 points to the 'View applied filters' button. A red box highlights the 'Filter List' panel at the bottom left, which displays 'Variant Type : SNP' and 'Consequence : missense_variant'. The main table on the right shows 213 variants of 74 variants, filtered by 'TEST'.

| Prediction | Tier | Warning | Low Conf... | Report | Gene | Transcript | Type | Consequence | NT Change |
|------------|------|---------|-------------|--------|--------|--------------|------|------------------|-----------|
| T1 | | | | | RUNX1 | NM_001754 | SNP | missense_variant | c.801G>A |
| T1 | | | | | RUNX1 | NM_001754 | SNP | missense_variant | c.146C>T |
| T2 | | | | | KRAS | NM_033360 | SNP | missense_variant | c.38G>A |
| T3 | | | | | DNMT1 | NM_001130823 | SNP | missense_variant | c.206G>A |
| T3 | | | | | CALR | NM_004343 | SNP | missense_variant | c.731A>G |
| T3 | | | | | NCKAP5 | NM_207363 | SNP | missense_variant | c.4117C>T |
| | | | | | CXCR4 | NM_001348056 | SNP | missense_variant | c.440A>G |
| | | | | | IDH1 | NM_001282386 | SNP | missense_variant | c.782C>T |
| | | | | | ASXL1 | NM_015338 | SNP | missense_variant | c.3745A>G |
| | | | | | EPC | NM_001126154 | SNP | missense_variant | c.1780C>C |

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

The screenshot shows the 'Variant Filter' dialog box from the NGeneAnalySys software. The 'Filter' dropdown menu is open, showing 'TEST' as the selected option. Below the dropdown are three tabs: 'Variant', 'Consequence', and 'Population Frequency'. The 'Variant' tab is active. Under the 'Variant' tab, there are several filter sections:

- Prediction:** Checkboxes for Tier1, Tier2, Tier3, and Tier4.
- Tier:** Checkboxes for Tier1, Tier2, Tier3, and Tier4.
- ClinVar:** Checkboxes for P(Pathogenic), LP(Likely Pathogenic), US(Uncertain Significance), LB(Likely Benign), and B(Benign).
- Variant Type:** Checkboxes for SNP, Ins(Insertion), Del(Deletion), MNP(Multiple-nucleotide polymorphism), and Complex(Multiple-nucleotides and an InDel).
- Chromosome:** A text input field containing 'chr7, chr9'.
- Gene:** A text input field containing 'FLT3,CEBPA'.
- dbSNP:** A checkbox labeled 'Has dbSNP ID'.
- COSMIC:** A checkbox labeled 'Has COSMIC ID'.
- COSMIC Occurrence:** A dropdown menu set to 'Tissue selection'.
- Allele Fraction:** Two input fields showing '30' and '50' with a percentage sign, separated by a tilde (~).

At the bottom right of the dialog box are two buttons: 'Cancel' and 'Save Filter', with 'Save Filter' being highlighted by a red box.

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The screenshot shows the NGeneAnalySys software interface with the following elements:

- Top Navigation:** DASHBOARD, RESULTS, SAMPLES (highlighted), MANAGER, Administrotor.
- Project Information:** HCT116-100, RUN: 2018-05-14 15:41:31, PANEL: HEMEaccuTest v1, DISEASE: Hematologic maligna..., SEQUENCER: Illumina MiSeqDx.
- Tab Selection:** OVERVIEW, VARIANT (highlighted with a yellow circle and number 1), REPORT.
- REPORTED VARIANT Section:** A table showing variants across different chromosomes (chr1, chr22) with details like Position, Transcript, NT Change, AA Change, and VAF. One variant is highlighted with a yellow circle and number 2.
- CONCLUSION Section:** A large empty text area highlighted with a yellow circle and number 3.
- TARGET GENES Section:** A grid of gene names (e.g., ABCB1, ARID1A, ATM, BCL11B, BLM, BCR, BRAF, BTG1, BTLA, CALN1, CARD11, CBL, CDKN1B, CDKN2A, CDKN2B, CHD1, CHD4, CREBBP, CRLF2, CSF1R). Some genes are highlighted with a yellow circle and number 4. Two specific genes, ARID1A and ATM, are highlighted with yellow circles and numbers 5 and 6 respectively.
- Action Buttons:** Gene Category dropdown (- N/A -), Save, Draft Report, Final Report.
- Version:** v 1.3 (Build Date 20180515).

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The screenshot shows the NGeneAnalySys software interface. At the top, there is a navigation bar with tabs: DASHBOARD, RESULTS, SAMPLES (highlighted with a yellow circle labeled 1), and MANAGER. Below the navigation bar, there is a sub-menu with tabs: Analysis Status, User Account, Panel, Report, News&Tips, System Logs, and Clinical Variant DB. A red box highlights the MANAGER tab. On the right side of the top bar, there is a user profile icon for 'Administrator'. The main content area displays a table of run information. The table has columns: Run ID, Run Name, Status, Login ID, User, Group, Date, Restart, and Delete. There are three rows of data:

| Run ID | Run Name | Status | Login ID | User | Group | Date | Restart | Delete |
|--------|-----------------|----------|----------|---------------|----------------|---------------------|---------|--------|
| 3 | BRCA_PLUS_FFPE | COMPLETE | admin | Administrator | Administrators | 2018-05-14 14:34:23 | | |
| 2 | BRCA_PLUS_BLOOD | COMPLETE | admin | Administrator | Administrators | 2018-05-14 14:34:07 | | |
| 1 | BRCA | COMPLETE | admin | Administrator | Administrators | 2018-05-14 14:33:54 | | |

At the bottom of the content area, there are navigation buttons: <<, ⏪, 1, ⏩, >>. To the right of these buttons, there is a 'Page' input field and a 'GO' button. A red box highlights the content area containing the table. In the bottom right corner of the window, there is a small text: 'v 1.3 (Build Date 20180515)'

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The screenshot shows the 'Analysis Status' section of the page. A red box highlights the 'Analysis Status' button at the top left. A yellow circle labeled '1' is over this button. Another yellow circle labeled '2' is over the table header. A third yellow circle labeled '3' is over the 'Refresh' button in the top right corner of the table area.

| Run ID | Run Name | Status | Login ID | User | Group | Date | Restart | Delete |
|--------|-----------------|----------|----------|---------------|----------------|---------------------|---------|--------|
| 3 | BRCA_PLUS_FFPE | COMPLETE | admin | Administrator | Administrators | 2018-05-14 14:34:23 | | |
| 2 | BRCA_PLUS_BLOOD | COMPLETE | admin | Administrator | Administrators | 2018-05-14 14:34:07 | | |
| 1 | BRCA | COMPLETE | admin | Administrator | Administrators | 2018-05-14 14:33:54 | | |

Below the table are navigation controls: '<<', '<', '1', '>', '>>', 'Page' (with a dropdown menu), and 'GO'.

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
The screenshot shows a user management interface with two main sections: 'Group' and 'User'. The 'User' section is highlighted by a large red box and contains three numbered callouts:

- 1**: A button labeled 'User Account' located at the top of the User panel.
- 2**: A section titled 'Group' containing a table with two rows. The first row has ID 2, Name 'Test', Deleted 'N', and Modify buttons. The second row has ID 1, Name 'Administra...', Deleted 'N', and Modify buttons.
- 3**: A section titled 'User' containing a table with two rows. The first row has Login ID 'Test', Name 'Test', User Type 'RESEARC...', Group 'Test', Created At '2018-06-04 04:38...', Updated At '-', Deleted 'N', and Modify buttons. The second row has Login ID 'admin', Name 'admin', User Type 'ADMIN', Group 'Administrators', Created At '2018-06-04 00:44...', Updated At '-', Deleted 'N', and Modify buttons.

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

| ID | Name | Deleted | Modify |
|----|----------------|---------|--------|
| 2 | Test | N | |
| 1 | Administrat... | N | |

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| ID | Name | Deleted | Modify |
|----|----------------|---------|---|
| 2 | Test | N |   |
| 1 | Administrat... | N |   |

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The screenshot shows a table titled "User" with one row of data. The columns are: Login ID, Name, User Type, Group, Created At, Updated At, Deleted, and Modify. The data in the first row is: admin, Adminis..., ADMIN, Administrators, 2018-05-17 08:26... - N. There are two red boxes: one around the "Modify" button in the top right corner labeled '1', and another around the entire row labeled '2'. A yellow circle with the number '3' is at the bottom right of the page.

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The screenshot shows the same User management interface as before, but the "Create" button at the bottom right is highlighted with a red box and a yellow circle with the number '3'.

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The screenshot shows the 'New User' dialog box from the NGeneAnalySys application. The dialog has a title bar 'NGeneAnalySys > New User'. The main area is titled 'New User'. It contains the following fields:

- Login ID (*) : Test
- Name (*) : Test
- Password (*) : [REDACTED]
- Confirm Password(*) : [REDACTED]
- User Type(*) : RESEARCHER
- User Group(*) : Test
- Organization : [REDACTED]
- Department : [REDACTED]
- Address : [REDACTED]
- Phone : [REDACTED]
- E-mail : [REDACTED]

At the bottom right are two buttons: 'Cancel' and 'Save', where 'Save' is highlighted in blue.

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
| User | | | | | | | |
|----------|------------|------------|----------------|---------------------|------------|---------|---|
| Login ID | Name | User Type | Group | Created At | Updated At | Deleted | Modify |
| Test | Test | RESEARC... | Test | 2018-06-04 04:38... | - | N |   |
| admin | Adminis... | ADMIN | Administrators | 2018-06-04 00:44... | - | N |   |

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 NGeneAnalySys > Edit User

Edit User

| | |
|--|------------|
| Login ID (*) | Test |
| Name (*) | Test |
| Password (*) | |
| Confirm Password (*) | |
| User Type (*) | RESEARCHER |
| User Group (*) | Test |
| Organization | |
| Department | |
| Address | |
| Phone | |
| E-mail | |
|   | |

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| User | | | | | | | |
|----------|------------|------------|----------------|---------------------|------------|---------|--------|
| Login ID | Name | User Type | Group | Created At | Updated At | Deleted | Modify |
| Test | Test | RESEARC... | Test | 2018-06-04 04:38... | - | N | |
| admin | Adminis... | ADMIN | Administrators | 2018-06-04 00:44... | - | N | |

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

The screenshot shows the 'Panel' tab selected in the top navigation bar. On the left, a red box highlights the 'Panel Update' configuration section for the 'SOLIDaccuTest v1' panel. This section includes fields for Panel name, Code, Target, Analysis Type, Library Type, Group, Disease, Default Sample Source, Default Disease, Report Template, Warning Read Depth(x), Essential Genes, Canonical Transcripts, Low Confidence, and Warning MAF(%). A blue checkmark icon with an 'Edit' label is present next to several of these fields. On the right, a table lists various panels with columns for Default, Panel name, Virtual Panel, Code, and Deleted status. The first panel in the list, 'SOLIDaccuTest v1', has its row highlighted with a light blue background. A red box surrounds the entire row of this panel, and another red box highlights the blue edit icons in the 'Virtual Panel' column for this row. Navigation buttons at the bottom indicate page 1 of 1.

| Default | Panel | Virtual Panel | Code | Deleted | |
|---------|------------------------|--|---------|---------|---|
| Y | SOLIDaccuTest v1 |   | - N/A - | 200 | N |
| Y | HEMEaccuTest v1 |   | - N/A - | 300 | N |
| Y | BRCAaccuTest v2 |   | - N/A - | 445 | N |
| Y | BRCAaccuTest PLUS v1 |   | - N/A - | 447 | N |
| Y | TruSight Tumor 170 DNA |   | - N/A - | 901 | N |
| Y | TruSight Tumor 170 RNA |   | - N/A - | 902 | N |

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The screenshot shows the 'Panel' tab selected in the top navigation bar. On the left, there is a 'Panel Update' sidebar with various configuration options. In the main area, a table lists several panels. The first panel, 'SOLIDaccuTest v1', has its delete button highlighted with a red box. The table columns include 'Default', 'Panel', 'Virtual Panel', 'Code', and 'Deleted'. The 'Deleted' column for the highlighted row shows 'N' (Not Deleted).

| Default | Panel | Virtual Panel | Code | Deleted |
|---------|------------------------|--|---------|---------|
| Y | SOLIDaccuTest v1 |   | - N/A - | 200 N |
| Y | HEMEaccuTest v1 |   | - N/A - | 300 N |
| Y | BRCAaccuTest v2 |   | - N/A - | 445 N |
| Y | BRCAaccuTest PLUS v1 |   | - N/A - | 447 N |
| Y | TruSight Tumor 170 DNA |   | - N/A - | 901 N |
| Y | TruSight Tumor 170 RNA |   | - N/A - | 902 N |

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


The screenshot shows the 'Report' tab selected in the top navigation bar. On the left, there's a detailed configuration panel for a new report. It includes fields for 'Report Name', 'Output Type', 'Parameter Setting' (with sub-fields for 'Name', 'Display Name', 'Type', and 'Combobox Item'), and sections for 'Parameters & Images' (with 'Parameters' and 'Images' dropdowns, and 'Choose Image File' and 'Choose Template File' buttons). At the bottom of this panel is a 'Save Parameter' button. On the right, a table lists existing reports with columns for 'Report ID', 'Report Name', 'Create Date', 'Update Date', 'Deleted', and 'Modify'. The first row in the table has a 'Modify' link highlighted with a red box and an arrow pointing to the 'Save Parameter' button in the left panel.

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

The screenshot shows the 'Report' tab selected in the top navigation bar. On the left, a red box highlights the 'Parameter Setting' section where report parameters like Name, Display Name, Type, and Combobox Item can be configured. On the right, a table lists existing reports with columns for Report ID, Report Name, Create Date, Update Date, Deleted status, and a 'Modify' button. The 'Modify' button for the first report ('testTemplate') is also highlighted with a red box.

| Report ID | Report Name | Create Date | Update Date | Deleted | Modify |
|-----------|--------------|-------------|-------------|---------|---|
| 1 | testTemplate | 2017-11-14 | 2018-02-19 | N |   |

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The screenshot shows a software interface for managing news tips. At the top, there are several tabs: Analysis Status, User Account, Panel, Report, News&Tips (which is highlighted with a red box and circled with a yellow circle labeled '1'), System Logs, and Clinical Variant DB. Below the tabs is a table with the following columns: ID, Title, Create Date, Update Date, Deleted, and Modify. The table has one row with ID 1, Title 'Test', Create Date '2018-05-17', Update Date '2018-05-17', Deleted 'N', and two icons for Modify. To the left of the table is a panel with a 'Title' input field (circled with a yellow circle labeled '2') and a 'Contents' area. At the bottom of this panel is a 'Save' button. To the right of the table is a red box containing a '+' button (circled with a yellow circle labeled '4')).

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The screenshot shows the 'News&Tips' tab selected in a navigation bar. On the left, there's a form for creating a new news item. It has fields for 'Title' (with a placeholder 'Title') and 'Contents'. A red box highlights this form area, and a red arrow points from the 'Save' button at the bottom right of the form to the '+ Add' button at the top right of the main content area. The main content area displays a table of existing news items:

| ID | Title | Create Date | Update Date | Deleted | Modify |
|----|-------|-------------|-------------|---------|--------|
| 1 | Test | 2018-05-17 | 2018-05-17 | N | |

At the bottom of the main content area, there are navigation icons: a left arrow, a page number '1', and a right arrow.

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The screenshot shows a user interface for managing news items. At the top, there is a navigation bar with several tabs: Analysis Status, User Account, Panel, Report, News&Tips (which is currently selected), System Logs, and Clinical Variant DB. Below the navigation bar, there is a table listing news items. The first news item in the list has its row highlighted with a red box. This row contains columns for ID (1), Title (Test), Create Date (2018-05-17), Update Date (2018-05-17), Deleted (N), and Modify. The 'Modify' column contains two icons: a blue edit icon and a red delete icon. A large red box encloses the entire left side of the screen, covering the navigation bar and the first news item's row. On the far right, there is a small blue '+' button. At the bottom of the screen, there is a footer with three circular icons and the number '1' in the center.

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The screenshot shows a software interface for managing news and tips. At the top, there are several tabs: Analysis Status, User Account, Panel, Report, News&Tips (which is selected and highlighted in blue), System Logs, and Clinical Variant DB. Below the tabs is a table with the following columns: ID, Title, Create Date, Update Date, Deleted, and Modify. There is one row in the table with the following data: ID 1, Title Test, Create Date 2018-05-17, Update Date 2018-05-17, Deleted N, and Modify (a button). The 'Modify' button is highlighted with a red box. To the left of the table, there is a vertical sidebar labeled 'Contents' and a 'Save' button at the bottom.

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The screenshot shows the 'System Logs' page. At the top, there is a navigation bar with several tabs: Analysis Status, User Account, Panel, Report, News&Tips, **System Logs**, and Clinical Variant DB. The 'System Logs' tab is highlighted with a red box and a yellow circle containing the number 1. Below the navigation bar, there is a search form with fields for Date, Log Type, Log Msg, User, and Group, followed by a 'Refresh' button. A red box surrounds the entire search form area, and a yellow circle containing the number 2 is placed over the 'Log Type' field. At the bottom of the page, there is a pagination control with buttons for <<, <, 1, >, >>, and a text input for 'Page' with a 'GO' button.

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
PATHOGENICITY CLASS

Pathogenic

- Nonsense mutations (stop_gained) with 1000 genomes and ESP6500 <= 0.01.
- Frameshift (frameshift_variant) with 1000 genomes and ESP6500 <= 0.01.
- Splice acceptor (splice_acceptor_variant) sites (-2 and -1) and/or splice donor (splice_donor_variant) sites (+1 and +2).
- Loss of natural initiation of transcription (start codon, start_lost) or termination (stop codon, stop_lost).
- In frame INDELs (inframe_insertion, inframe_deletion) of more than 3 nucleotides.
- Mutation reported as pathogenic in BIC or ClinVar.

Likely Pathogenic

- Missense (missense_variant) mutation with 1000 genomes and ESP6500 < 0.01.
- In frame INDELs of 3 nucleotides.
- Missense or synonymous with at least 2 values of predictive score (SIFT, PolyPhen2, MutationTaster) > 0.9.

Likely Benign

- None.

Benign

- Missense, synonymous or intronic variant with 1000 genomes or ESP6500 > 0.05 with 2 predictive.

Uncertain Significance

- Mutation no classified in pathogenic, likely pathogenic and benign.

TIER CLASS

Tier I

- Variants with strong clinical significance.
- Level A: biomarkers that predict response or resistance to US FDA-approved therapies for a specific type of tumor or have been included in professional guidelines as therapeutic, diagnostic, and/or prognostic biomarkers for specific types of tumors.
- Level B: biomarkers that predict response or resistance to a therapy based on well-powered studies with consensus from experts in the field, or have diagnostic and/or prognostic significance of certain diseases based on well powered studies with expert consensus.
- Variants with strong clinical significance.

Tier II

- Variants with potential clinical significance.
- Level C: biomarkers that predict response or resistance to therapies approved by FDA or professional societies for a different tumor type (ie, off-label use of a drug), serve as inclusion criteria for clinical trials, or have diagnostic and/or prognostic significance based on the results of multiple small studies.
- Level D: biomarkers that show plausible therapeutic significance based on preclinical studies, or may assist disease diagnosis and/or prognosis themselves or along with other biomarkers based on small studies or multiple case reports with no consensus.

Tier III

- Variants of unknown clinical significance.
- Not observed at a significant allele frequency in the general or specific subpopulation databases, or pan-cancer or tumor-specific variant databases.
- No convincing published evidence of cancer association.

Tier IV

- Variants deemed benign or likely benign.
- Observed at significant allele frequency in the general or specific subpopulation databases.
- No existing published evidence of cancer association.

PATHOGENICITY CLASS - ClinVar

Clinical significance on ClinVar submissions (SCV)

| Clinical Significance | | Guidance for Use in ClinVar SCV Records |
|---|-------|---|
| | Value | |
| Benign | | As recommended by ACMG/AMP for variants interpreted for Mendelian disorders. |
| Likely benign | | As recommended by ACMG/AMP for variants interpreted for Mendelian disorders. |
| Uncertain significance | | As recommended by ACMG/AMP for variants interpreted for Mendelian disorders. |
| Likely pathogenic | | As recommended by ACMG/AMP for variants interpreted for Mendelian disorders. |
| Pathogenic | | As recommended by ACMG/AMP for variants interpreted for Mendelian disorders. |
| drug response | | A general term for a variant that affects a drug response, not a disease. We anticipate adding more specific drug response terms based on a recommendation by CPIC. |
| association | | For variants identified in a GWAS study and further interpreted for their clinical significance. |
| risk factor | | For variants that are interpreted not to cause a disorder but to increase the risk. |
| protective | | For variants that decrease the risk of a disorder, including infections. |
| Affects | | For variants that cause a non-disease phenotype, such as lactose intolerance. |
| conflicting data from submitters | | Only for submissions from a consortium, where groups within the consortium have conflicting interpretations of a variant but provide a single submission to ClinVar. |
| other | | If ClinVar does not have the appropriate term for your submission, we ask that you submit "other" as clinical significance and contact us to discuss if there are other terms we should add. |
| not provided | | For submissions without an interpretation of clinical significance. The primary goal of ClinVar is to archive reports of clinical significance of variants. Therefore, submissions with a clinical significance of "not provided" should be limited to: <ul style="list-style-type: none"> • "literature only" submissions that report a publication about the variant, without interpreting the clinical significance. • "research" submissions that provide functional significance (e.g. |

| | |
|--|---|
| | <p>undetectable protein level) but no interpretation of clinical significance.</p> <ul style="list-style-type: none"> “clinical testing” or “phenotyping only” submissions from clinics or physicians that provide additional information about individuals with the variant, such as observed phenotypes, but do not interpret the clinical significance. |
|--|---|

(<https://www.ncbi.nlm.nih.gov/clinvar/docs/clinsig/>)

Clinical significance on ClinVar aggregate records (RCV)

| Combination of Values from Different Submitters | Reported as |
|---|---|
| Pathogenic and Likely pathogenic | Pathogenic/Likely Pathogenic |
| (Pathogenic or Likely pathogenic or Benign or Likely benign) AND Uncertain significance | Conflicting interpretations of pathogenicity |
| (Pathogenic or Likely pathogenic) AND (Benign or Likely benign) | Conflicting interpretations of pathogenicity |
| Benign and Likely benign | Benign/Likely benign |
| Any single ACMG value AND any non-ACMG value, e.g. Uncertain significance and risk factor | Uncertain significance, risk factor |
| Conflicting ACMG values AND any non-ACMG value, e.g. Pathogenic and Uncertain significance and risk factor | Conflicting interpretations of pathogenicity, risk factor |
| No ACMG value AND multiple non-ACMG values, e.g. Drug response and risk factor | drug response, risk factor |
| A consortium does not agree on which ACMG assessment is valid | conflicting data from submitters |

(<https://www.ncbi.nlm.nih.gov/clinvar/docs/clinsig/>)

PATHOGENICITY CLASS - BIC

| Class | Clinical Significance Value | Description |
|----------------|--|---|
| Class 1 | Not pathogenic low clinical significance | There is significant evidence against this variant being a dominant high-risk pathogenic mutation. |
| Class 2 | Likely Notice pathogenic little clinical significance | There is evidence against this variant being a dominant high-risk pathogenic mutation. |
| Class 3 | Uncertain | There is insufficient evidence to place this variant in Class 1,2,4 or 5. |
| Class 4 | Likely Pathogenic | There is evidence that this variant is a dominant high-risk pathogenic mutation. |
| Class 5 | Pathogenic | There is significant evidence to suggest that this variant is a dominant high-risk pathogenic mutation. |

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| 번호 | 필드명 | 패널 | Default | 설명 |
|-------------------------------|--------------------|--------------|---------|---|
| 사용자 설정 항목 | | | | |
| 1 | Prediction | All | Y | Pathogenic class using ACMG ISV or AMP/ASCO.CAP guideline. |
| 2 | Pathogenicity | Germline DNA | Y | User defined Pathogenic class. |
| 3 | Tier | Somatic DNA | Y | User defined Tier class. |
| 4 | Warning | All | Y | Warning variant. |
| 5 | Low Confidence | Somatic DNA | Y | Low confidence variant |
| 6 | Report | All | Y | Reported variant. |
| 변이 및 유전자 annotation 정보 | | | | |
| 7 | Gene | All | Y | Gene Symbol (HUGO Gene Nomenclature). |
| 8 | Strand | All | N | Whether the gene is located on the plus or minus strand. |
| 9 | Transcript | All | Y | Transcript ID. |
| 10 | Type | All | Y | Variant type SNP: Single-Nucleotide Polymorphism, INS: Insertion, Del: Deletion, MNP: Multiple-nucleotide polymorphism, Complex: Multiple-nucleotide and an Indel. |
| 11 | Type Extension | Somatic DNA | Y | |
| 12 | Consequence | All | Y | Sequence change result by variant. |
| 13 | NT Change | All | Y | The variant coordinates relative to the coding DNA |
| 14 | NT Change (BIC) | Germline DNA | Y | The variant coordinates relative to the coding DNA (BRCA format-BRCA1: A of start codon ATG is #120, BRCA2: A of start codon ATG is #229). |
| 15 | AA Change (Single) | All | Y | Change amino acid (single letter). |
| 16 | AA Change | All | Y | Change amino acid. (three letter) |
| 17 | dbSNP ID | All | Y | dbSNP RS ID. |
| 18 | dbSNP Common ID | All | N | dbSNP RS Common ID. |
| 19 | Exon | All | Y | Number of exon in which the variant is located exon/total exon. |
| 20 | Exon (BIC) | Germline DNA | Y | Exon number; Reference sequences: BRCA1 GenBank U14680; BRCA2 GenBank U43746. |

| 변이 위치 정보 | | | | |
|-------------------|-----------------------|--------------|---|--|
| 21 | Chr | All | Y | The chromosome numbers. |
| 22 | Start Position | All | Y | The variant coordinate on the reference genome. |
| 23 | Ref | All | Y | Reference allele. |
| 24 | Alt | All | Y | Alternative allele. |
| 25 | Zigosity | Germline DNA | N | The words homozygous, heterozygous, and hemizygous are used to describe the genotype of a diploid organism at a single locus on the DNA. |
| Allele 정보 | | | | |
| 26 | Fraction | All | Y | The proportion of the reads supporting the variant count. |
| 27 | Depth | All | Y | Total coverage of the variant. |
| 28 | Ref Count | All | Y | Number of reads supporting the reference allele. |
| 29 | Alt Count | All | Y | Number of reads supporting the alternative allele. |
| ClinVar 데이터베이스 정보 | | | | |
| 30 | ClinVar Accession | All | Y | ClinVar RCV accessions. |
| 31 | ClinVar Class | All | Y | Clinical classification. |
| 32 | ClinVar Disease | All | Y | ClinVar disease name. |
| 33 | ClinVar Trait OMIM | All | Y | |
| Population 빈도 | | | | |
| 34 | 1KGP All | All | Y | 1000 Genomes Project(1KGP)'s All population frequency. |
| 35 | 1KGP African | All | N | 1000 Genomes Project(1KGP)'s African population frequency. |
| 36 | 1KGP American | All | N | 1000 Genomes Project(1KGP)'s American population frequency. |
| 37 | 1KGP East Asian | All | Y | 1000 Genomes Project(1KGP)'s East Asian population frequency. |
| 38 | 1KGP European | All | N | 1000 Genomes Project(1KGP)'s European population frequency. |
| 39 | 1KGP South Asian | All | N | 1000 Genomes Project(1KGP)'s South Asian population frequency. |
| 40 | ESP All | All | Y | NHLBI Exome Sequencing Project(ESP)'s All population frequency. |
| 41 | ESP African American | All | N | NHLBI Exome Sequencing Project(ESP)'s African American population frequency. |
| 42 | ESP European American | All | N | NHLBI Exome Sequencing Project(ESP)'s European American population frequency. |
| 43 | ExAC | All | Y | Exome Aggregation Consortium(ExAC) population frequency. |
| 44 | gnomAD All | Somatic DNA | Y | Genome Aggregation Database(gnomAD)'s All population frequency. |

| | | | | |
|-----------|---------------------------------|-------------|---|--|
| 45 | gnomAD Admixed American | Somatic DNA | N | Genome Aggregation Database(gnomAD)'s Admixed American population frequency. |
| 46 | gnomAD African African American | Somatic DNA | N | Genome Aggregation Database(gnomAD)'s African African American population frequency. |
| 47 | gnomAD East Asian | Somatic DNA | Y | Genome Aggregation Database(gnomAD)'s East Asian population frequency. |
| 48 | gnomAD Finnish | Somatic DNA | N | Genome Aggregation Database(gnomAD)'s Finnish population frequency. |
| 49 | gnomAD Non Finnish European | Somatic DNA | N | Genome Aggregation Database(gnomAD)'s Non Finnish European population frequency. |
| 50 | gnomAD Others | Somatic DNA | N | Genome Aggregation Database(gnomAD)'s Others population frequency. |
| 51 | gnomAD South Asian | Somatic DNA | N | Genome Aggregation Database(gnomAD)'s South Asian population frequency. |
| 52 | KRGDB | All | Y | Korean Reference Genome Database(KRGDB) population frequency. |
| 53 | KoEXID | All | N | Korean Exome Information Database(KoEXID) population frequency. |

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| | | | | |
|-----------|-----------------|--------------|---|---|
| 54 | BIC Category | Germline DNA | N | Mutation type in BIC. 3UTR: 3'UTR. 5TR: 5' UTR. F: Frameshift. IFI: In frame insertion. IFD: In frame deletion. IVS: Intervening sequence. M: Missense. N: Nonsense. P: Polymorphism. S: Splice. Syn: Synonymous. UV: Unclassified variant. |
| 55 | BIC Class | Germline DNA | Y | Clinical classification using 5 classes in BIC. 0: pending. 1: Not pathogenic/low clinical significance. 2: Likely not pathogenic/little clinical significance. 3: Uncertain. 4: Likely pathogenic. 5: Pathogenic. |
| 56 | BIC Designation | Germline DNA | Y | Designation of described mutation according to BIC nomenclature guidelines. |
| 57 | BIC Importance | Germline DNA | N | Clinical importance in BIC. |

| | | | | |
|-----------|--------------------------|--------------|---|---|
| 58 | BIC NT | Germline DNA | Y | Nucleotide # in the transcript (cDNA) at which mutation occurs; Reference sequences: BRCA1 GenBank U14680; BRCA2 GenBank U43746. |
| 59 | KOHBRA Frequency | Germline DNA | Y | High risk breast cancer frequency in Korea (KOHBRA). |
| 60 | KOHBRA Patient | Germline DNA | Y | High risk breast cancer patient in Korea (KOHBRA). |
| 61 | Be BIC Category | Germline DNA | N | Variants are classified as either germline or somatic, depending on how they are acquired. Germline variants are genetic changes that we inherit from our parents. Somatic variants are DNA changes that we acquire over our lifetime, often through exposure to pollutants, toxins, radiation and other carcinogens. |
| 62 | Be BIC Ethnic | Germline DNA | N | Method used to collect the data that supports the assertion of clinical significance. Allowed values: case-control, clinical testing, literature only, reference population, research. |
| 63 | Be BIC Nationality | Germline DNA | N | Mutation type in BIC. |
| 64 | Be BIC Pathogenicity | Germline DNA | N | Nationality of the patient associated with the variant. From BIC. |
| 65 | Be ClinVar Method | Germline DNA | N | Ethnic of BIC. |
| 66 | Be ClinVar Origin | Germline DNA | N | BIC pathogenicity. |
| 67 | Be ClinVar Pathogenicity | Germline DNA | N | Reference transcript. |
| 68 | Be ClinVar Update | Germline DNA | N | Date the variant was last updated in ClinVar. |
| 69 | Be ENIGMA Condition | Germline DNA | N | Nucleotide change. |
| 70 | Be ENIGMA Pathogenicity | Germline DNA | N | Reference gene. |
| 71 | Be ENIGMA Update | Germline DNA | N | Human-readable condition, describing the biological impact of the variation. |
| 72 | Be Gene | Germline DNA | N | The date on which the clinical significance of the variant was last evaluated by ENIGMA. |
| 73 | Be NT | Germline DNA | N | IRAC Classification from exLOVD. |
| 74 | Be Transcript | Germline DNA | N | ClinVar pathogenicity. |
| 75 | ENIGMA | Germline DNA | N | The Clinical Significance columns indicates whether expert curators have determined if the variant is |

| | | | | |
|----|----------------------|-------------|---|--|
| | | | | pathogenic or benign. What do these classifications mean? Pathogenic: variants confer an increased risk of disease. Likely pathogenic: variants have good evidence to support an association with disease risk. Likely benign: variants have good evidence to support no association with disease risk. Benign: variants are not associated with any markedly increased risk of disease. Variants of uncertain significance (VUS): are those for which the evidence of disease risk is not clear yet, sometimes because there is not yet enough evidence to classify them as either pathogenic or benign. |
| 76 | COSMIC ID | Somatic DNA | N | COSMIC ID. |
| 77 | COSMIC Occurrence | Somatic DNA | N | Tissue classification in COSMIC. |
| 78 | COSMIC Count | Somatic DNA | N | Sample count detected COSMIC variant. |
| 79 | OncoKB Action Cancer | Somatic DNA | N | OncoKB actionable gene : cancer information |
| 80 | OncoKB Action Drugs | Somatic DNA | N | OncoKB actionable gene : drug information |
| 81 | OncoKB Action Level | Somatic DNA | N | OncoKB actionable gene : evidence level; R1:Standard of care biomarker predictive of resistance to an FDA-approved drug in this indication / R2:Compelling clinical evidence supports the biomarker as being predictive of resistance to a drug, but neither biomarker nor drug are standard care / R3:Compelling biological evidence supports the biomarker as being predictive of resistance to a drug, but neither biomarker nor drug are standard care / 1:FDA-approved biomarker and drug in this indication / 2A:Standard of care biomarker predictive of response to an FDA-approved drug in this indication / 2B:Standard of care biomarker predictive of response to an FDA-approved drug in another indication but not standard of care for this indication / 3A:Compelling clinical evidence supports the biomarker as being predictive of response to a drug in this indication but neither biomarker and drug are standard of care / 3B:Compelling clinical evidence supports the biomarker as being predictive of response to a drug in another indication but neither |

| | | | | |
|-----------|---------------------|-------------|---|---|
| | | | | biomarker and drug are standard of care / 4:Compelling biological evidence supports the biomarker as being predictive of response to a drug but neither biomarker and drug are standard of care |
| 82 | OncoKB Action PMID | Somatic DNA | N | OncoKB actionable gene : PMIDs |
| 83 | OncoKB Alteration | Somatic DNA | N | OncoKB annotated gene : Mutation effect |
| 84 | OncoKB HGVS.p | Somatic DNA | N | OncoKB annotated gene : HGVS.p |
| 85 | OncoKB Oncogenicity | Somatic DNA | N | OncoKB annotated gene : Oncogenicity |
| 86 | OncoKB PMIDs | Somatic DNA | N | OncoKB annotated gene : PMIDs |
| 87 | OncoKB TX | Somatic DNA | N | OncoKB annotated gene : Transcript number |

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- 1) Lynch HT, Snyder CL, Lynch JF, Riley BD, Rubinstein WS. Hereditary breast-ovarian cancer at the bedside: role of the medical oncologist. *Journal of clinical oncology*. 2003;21(4):740-53.
- 2) Robson M, Offit K. Clinical practice. Management of an inherited predisposition to breast cancer. *The New England journal of medicine*. 2007;357(2):154-62.
- 3) Pal T, Permuth-Wey J, Betts JA, et al. BRCA1 and BRCA2 mutations account for a large proportion of ovarian carcinoma cases. *Cancer*. 2005;104(12):2807-16.
- 4) Kang E, Seong MW, Park SK, et al. The prevalence and spectrum of BRCA1 and BRCA2 mutations in Korean population: recent update of the Korean Hereditary Breast Cancer (KOHBRA) study. *Breast Cancer Res Treat*. 2015;151:157-68.
- 5) Moyer VA. Risk assessment, genetic counseling, and genetic testing for BRCA-related cancer in women: U.S. Preventive Services Task Force recommendation statement. *Annals of internal medicine*. 2014;160(4):271-81.
- 6) Gadzicki D, Evans DG, Harris H, et al. Genetic testing for familial/hereditary breast cancer-comparison of guidelines and recommendations from the UK, France, the Netherlands and Germany. *Journal of community genetics*. 2011;2(2):53-69.
- 7) American Society of Clinical Oncology policy statement update: genetic testing for cancer susceptibility. *Journal of clinical oncology*. 2003;21(12):2397-406.
- 8) Rechards S, Aziz N, Bale S, et al. Standards and guidelines for the interpretation of sequence variants: a joint consensus recommendation of the American College of Medical Genetics and Genomics and the Association for Molecular Pathology. *Genetics in Medicine*. 2015;17(5):405-24.
- 9) CPIC Term Standardization for Clinical Pharmacogenetic Test Results Project.
<https://www.pharmgkb.org/page/cpicTermProject>
- 10) Representation of clinical significance in ClinVar and other variation resources at NCBI.
<https://www.ncbi.nlm.nih.gov/clinvar/docs/clinsig/>

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*Copyright © 2016 by NGeneBio Co., Ltd. All Rights Reserved.