NGeneAnalySys™

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<체외진단용 의료기기>

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
The screenshot shows the main interface of NGeneAnalySys. At the top, there are tabs for 'HOME', 'Past Results', and 'Patient DB'. On the right, there's a red 'Upload FASTQ Files' button. The main area is titled 'Running and Recent Samples' with a sub-section 'Running and Recent Samples (Running : 0 / Error : 0)'. It lists four recent samples:

Sample	Run	Job Status	Panel	QC
B2_S2_L001	2017-06-21 20:11:36 Submitted: 2017-06-21 20:11:36 Started: 2017-06-21 20:13:24	Specimen: BLOOD REPORTING	UPLOAD DATA → RUNNING ANALYSIS → COMPLETE NGeneBio BRCAccuTest® Illumina MiSeq	PASS
B1_S1_L001	2017-06-21 20:10:45 Submitted: 2017-06-21 20:10:45 Started: 2017-06-21 20:13:26	Specimen: BLOOD PIPELINE FINISHED	UPLOAD DATA → RUNNING ANALYSIS → COMPLETE NGeneBio BRCAccuTest® Illumina MiSeq	PASS
NA12878_S01_L001	2017-06-14 14:16:54 Submitted: 2017-06-14 14:16:54 Started: 2017-06-14 14:16:54	Specimen: BLOOD PIPELINE FINISHED	UPLOAD DATA → RUNNING ANALYSIS → COMPLETE NGeneBio BRCAccuTest® Illumina MiSeq	PASS
NA12878_S01_L001	2017-06-14 13:16:45 Submitted: 2017-06-14 13:16:45 Started: 2017-06-14 14:13:22	Specimen: BLOOD PIPELINE FINISHED	UPLOAD DATA → RUNNING ANALYSIS → COMPLETE NGeneBio BRCAccuTest® Illumina MiSeq	PASS


At the bottom, there's a 'Progress Task' bar and the text 'v1.1 (Build Date 20170530)'.

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
The screenshot shows the NGeneAnalySys software interface. At the top, there is a header bar with the NGeneBio logo and navigation links for HOME, Past Results, and Patient DB. Below the header is a search bar. The main content area is titled "Run Analysis" and displays a message about the BRCAAccuTest panel. A large red box highlights the "Upload FASTQ Files" button (labeled 2). Another red box highlights the "Running and Recent Samples" table (labeled 1). A third red box highlights the "Progress Task" section at the bottom (labeled 3).

Sample	Run	Job Status	Panel	QC
B2_S2_L001 Submitted : 2017-06-21 20:11:36 Started : 2017-06-21 20:13:24	2017-06-21-20-11 Specimen : BLOOD	UPLOAD DATA > RUNNING ANALYSIS > COMPLETE REPORTING	NGeneBio BRCAAccuTest® Illumina MiSeq	PASS
B1_S1_L001 Submitted : 2017-06-21 20:10:45 Started : 2017-06-21 20:13:26	2017-06-21-20-10 Specimen : BLOOD	UPLOAD DATA > RUNNING ANALYSIS > COMPLETE PIPELINE FINISHED	NGeneBio BRCAAccuTest® Illumina MiSeq	PASS
NA12878_S01_L001 Submitted : 2017-06-14 14:16:11 Started : 2017-06-14 14:16:54	2017-06-14-14-16 Specimen : BLOOD	UPLOAD DATA > RUNNING ANALYSIS > COMPLETE PIPELINE FINISHED	NGeneBio BRCAAccuTest® Illumina MiSeq	PASS
NA12878_S01_L001 Submitted : 2017-06-14 13:16:45 Started : 2017-06-14 14:13:22	2017-06-14-13-16 Specimen : BLOOD	UPLOAD DATA > RUNNING ANALYSIS > COMPLETE PIPELINE FINISHED	NGeneBio BRCAAccuTest® Illumina MiSeq	PASS

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The screenshot shows the 'New Analysis Request' dialog box. The top section contains form fields for basic analysis parameters, each with a question mark icon and a dropdown or input field. A green circle labeled '1' is positioned to the left of this section. Below the form is a table with four columns: MID, Sample, Experiment Type, and Panel. The table is currently empty and displays the message 'No content in table'. A green circle labeled '2' is positioned to the left of the table. At the bottom right of the dialog are two buttons: 'Cancel' (gray) and 'Next' (red), with a green circle labeled '3' positioned above the 'Next' button.

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
MID	Sample	Experiment Type	Panel
S19	NA14094-1	germline	NGeneBio BRCAccuTest®
S3	NA13707-1	germline	NGeneBio BRCAccuTest®

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Running and Recent Samples (Running : 0 / Error : 0)		
Sample	Run	Job Status
B2_S2_L001 Submitted : 2017-06-21 20:11:36 Started : 2017-06-21 20:13:24	2017-06-21-20-11 Specimen : BLOOD	UPLOAD DATA > RUNNING ANALYSIS > COMPLETE REPORTING
B1_S1_L001 Submitted : 2017-06-21 20:10:45 Started : 2017-06-21 20:13:26	2017-06-21-20-10 Specimen : BLOOD	UPLOAD DATA > RUNNING ANALYSIS > COMPLETE PIPELINE FINISHED

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UPLOAD DATA > RUNNING ANALYSIS > COMPLETE In preparation (yellow)
Uploading waiting

UPLOAD DATA > RUNNING ANALYSIS > COMPLETE In progress (green)

UPLOAD DATA > RUNNING ANALYSIS > COMPLETE Complete of status (blue)
[STATUS] FASTQ file validation

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The screenshot shows the NGeneBio software interface. At the top, there is a navigation bar with tabs: HOME, Past Results (highlighted with a red box and circled with number 1), and Patient DB. Below the navigation bar is a search bar labeled "Sample Filter". The main content area displays a table titled "Completed Samples (4 Samples)". The table has columns for Sample, Run, Panel, and Result Overview. Each row contains sample details, run information, panel name (e.g., "NGeneBio BRCAaccuTest® Illumina MiSeq"), and a result summary (e.g., GENES: 2, DEPTH MIN: 189, DEPTH MAX: 564). A "Detail" button is present in each row. The bottom right corner of the interface shows the version "v 1.1 (Build Date 201705303)".

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The screenshot shows the "Sample Filter" search interface. It includes fields for "Panel", "Specimen", "Run", "Sample", "Submitted", and "Status". A "Search" button is located on the right side of the filter bar (circled with number 3). A "Reset" button is also present (circled with number 2).

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분류	아이콘	설명
Pathogenic	P	Causative for clinical phenotype.
Likely Pathogenic	LP	Probably causative for clinical phenotype.
Uncertain Significance	US	Unknown effect on clinical phenotype.
Likely Benign	LB	Probably not causative for a clinical phenotype.
Benign	B	Not causative for a clinical phenotype.

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Prediction	Pathogenicity	Warn	Report	False	Type	Cod.Cons	Gene	Strand	Transcript	NT change
B					R		BRCA2	+	NM_000059.3	c.1114A>C
B					R		BRCA2	+	NM_000059.3	c.3396A>G
B					snp	synonymous	BRCA2	+	NM_000059.3	c.4563A>G
B					snp	synonymous	BRCA2	+	NM_000059.3	c.6513G>C

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

Visualization for variant location

Provide extra window to see variant information

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
Variant List - sorted by:					
Prediction	Pathogenicity	Warn	Report	False	
B					
B					
P	LP	▲	R	F	
B					
B					

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The screenshot shows a table with columns: Type, Filename, Size, Created, and Download. A red box highlights the first row. Green circles numbered 1 through 6 point to specific elements: 1 points to the 'RAW Data' tab; 2 points to the 'Type' column header; 3 points to the 'Filename' column header; 4 points to the 'Size' column header; 5 points to the 'Created' column header; 6 points to the 'Download' button for the QC_Report row.

Type	Filename	Size	Created	Download
FASTQ	s5-2_S5_L001_R2_001.fastq.gz	2.7 MB	2017-02-02 19:32:47	Download
FASTQ	s5-2_S5_L001_R1_001.fastq.gz	2.6 MB	2017-02-02 19:32:47	Download
BAM	s5-2_S5_L001_final.bam	1.4 MB	2017-02-02 19:32:47	Download
BAM	s5-2_S5_L001_final.bai	36 kB	2017-02-02 19:32:47	Download
VCF	s5-2_S5_L001.vcf	3 kB	2017-02-02 19:32:47	Download
QC_Report	s5-2_S5_L001.pdf	2.1 MB	2017-02-02 19:32:47	Download

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The screenshot shows the 'Reports' tab of a clinical reporting software. The interface is divided into several sections:

- Test Information:** Includes a dropdown menu for 'Recorded Inspection Item' (2), a 'Manage No' input field, and an 'Inspection Item' input field.
- Client Information:** Includes a dropdown menu for 'Recorded Client Info' (3), fields for 'Organization', 'Client (Doctor)', and 'Contact'.
- Patient Information:** Includes fields for 'Patient ID', 'Name' (first name and last name), 'Birthday' (yyyy-MM-dd format), 'Gender' (All), and buttons for 'Add To Patient DB' and 'Search' (4).
- Reported Variant List:** A table header with columns: Prediction, Pathogenicity, Warn, False, Type, Cod.Cons, Gene, Transcript, NT change, AA.change, NT.change(BIC), Chr, Ref, Alt, Zygosity, g.pos, Exon, Exon(BIC), ref.num, alt.n. Below the table is a message: 'No data in table' (5).
- Conclusions:** An empty text area (6).

At the bottom right are three buttons: 'Save' (7), 'Create Draft Report', and 'Create Final Report'.

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The screenshot shows the 'Edit' dialog box with two main sections:

- Change Password** (Section 1):
 - Current Password
 - New Password
 - New Password Confirm
 - Hint: Password is English, numbers, please enter at least 8 or more characters with a combination of special characters.
- User Information** (Section 2):
 - Organization
 - Department
 - Address
 - Phone
 - E-mail

At the bottom right are the **Save** and **Close** buttons.

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
The screenshot shows a software window titled "NGeneAnalySys > Public Databases". The window contains a table with columns: Database, Version, and Source. The table lists various genomic databases and their sources:

Database	Version	Source
1000 Genomes	Phase3 (May, 2013)	ftp://ftp.1000genomes.ebi.ac.uk/vol1/ftp/release/20130502/
ExAC	Release 1.0 (Feb27, 2017)	ftp://ftp.broadinstitute.org/pub/ExAC_release/release1/ExAC.r1.sites.vep.vcf.gz
dbSNP	dbSNP 150 (Feb3,2017)	ftp://ftp.ncbi.nih.gov/snp/organisms/human_9606_b150_GRCh37p13/VCF/00-All.vcf.gz
ESP6500	ESP6500SI-V2	http://evs.gs.washington.edu/evs_bulk_data/ESP6500SI-V2-SSA137.GRCh38-liftOver.snps_indels.vcf.tar.gz
Reference genome	GRCh37/hg19	
Clinvar	Sep, 2017	ftp://ftp.ncbi.nlm.nih.gov/pub/clinvar/xml/ClinVarFullRelease_2017-09.xml.gz
BIC	Oct, 2017	https://research.nhgri.nih.gov/bic/
dbNSFP	v3.5 (Aug6,2017)	https://sites.google.com/site/jpopgen/dbNSFP
gnomAD	v2.0.1 (Feb27, 2017)	http://gnomad.broadinstitute.org/downloads
KOHBRA	Breast Cancer Res Treat(2015) 151:157-168	Breast Cancer Res Treat(2015) 151:157-168
KOEX	Experimental & Molecular Medicine (2017) 49, e356	Experimental & Molecular Medicine (2017) 49, e356
ARUP	April, 2017	http://arup.utah.edu/database/BRCA/
KRGDB	krgdb_622(2014), krgdb_1100(2016)	http://152.99.75.168/KRGDB/menuPages/download.jsp
dbscSNV	v1.1 (Apr12, 2015)	ftp://dbnsfp:dbnsfp@dbnsfp.softgenetics.com/dbscSNV.zip

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The screenshot shows the 'User Management' section of the NGeneAnalySys Manager. The left sidebar has a 'User Management' menu item highlighted with a red box and a green circle labeled '1'. The main content area displays a table titled 'User list' with columns: ID, LOGIN_ID, NAME, USER_TYPE, GROUP_ID, GROUP_NAME, CREATED_AT, UPDATED_AT, DELETED_AT, DELETED, and USER MODIFY. The table contains three rows of data. At the top right of the content area, there are buttons for 'Create User' and 'Delete User', and a 'Logout' button.

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The screenshot shows the application's main menu on the left and the User Management page on the right. The User Management page has two main sections: 'Group list' and 'User list'. Both sections include search and filter fields. The 'User list' section displays a table of user records with columns: ID, LOGIN_ID, NAME, USER_TYPE, GROUP_ID, GROUP_NAME, CREATED_AT, UPDATED_AT, DELETED_AT, and USER MODIFY. Three user entries are shown: test (RESEARCHER), controller (ADMIN), and admin (ADMIN).

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The screenshot shows a 'Group list' interface. At the top right are two buttons: '+ Create Group' and '- Delete Group'. A red box surrounds these buttons. Below them is a search bar with 'Group Name:' and a 'Search' button; a red box surrounds this area. The main content is a table with columns: ID, GROUP_NAME, DELETED, and GROUP MODIFY. It contains one row with ID 1, GROUP_NAME 'administrators', and DELETED 'N'. A red box surrounds the entire table. Three green circles with numbers 1, 2, and 3 are overlaid on the interface: 1 is on the '+ Create Group' button, 2 is on the search bar, and 3 is on the table body.

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The screenshot shows the same 'Group list' interface as before, but the '+ Create Group' button is highlighted with a red box. Below it, a 'Group Create' dialog box is open. The dialog has a title 'Group Create' and a single input field labeled 'Group Name' with the placeholder 'Group Name'. At the bottom are 'Close' and 'Save' buttons.

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A close-up screenshot of the 'Group Create' dialog box. It shows the 'Group Name' input field containing 'Group Name'. At the bottom are two buttons: 'Close' on the left and 'Save' on the right.

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Group list			
	 Create Group	 Delete Group	
ID	GROUP_NAME	DELETED	GROUP MODIFY
<input type="checkbox"/>	1 administrators	N	

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Group list				
		+ Create Group	- Delete Group	
Group Name :				Search
ID	GROUP_NAME	DELETED	GROUP MODIFY	
<input checked="" type="checkbox"/>	2 Hospital	N		
<input type="checkbox"/>	1 administrators	N		
				1

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Group list				
		+ Create Group	- Delete Group	
Group Name :		Hospital	Search	
ID	GROUP_NAME	DELETED	GROUP MODIFY	
<input type="checkbox"/>	2 Hospitals	N		
				1

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User list

ID	LOGIN_ID	NAME	USER_TYPE	GROUP_ID	GROUP_NAME	CREATED_AT	UPDATED_AT	DELETED_AT	DELETED	USER MODIFY
3	test	test	RESEARCHER	1	administrators	2017-04-28T01:53:28.396955Z	-	-	N	
2	controller	controller	ADMIN	1	administrators	2016-05-19T00:55:30.748446Z	-	-	N	
1	admin	admin	ADMIN	1	administrators	2016-05-19T00:55:30.748446Z	-	-	N	

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+ Create User

Type : ALL	Group Name : ALL	Search			
NAME	CREATED_AT	UPDATED_AT	DELETED_AT	DELETED	USER MODIFY
ors	2017-04-28T01:53:28.396955Z	-	-	N	
ors	2016-05-19T00:55:30.748446Z	-	-	N	
ors	2016-05-19T00:55:30.748446Z	-	-	N	

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User Create

Login ID	Login ID
Name	Name
Password	Password
User Type	ADMIN
User Group	Hospitals

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User list										
Login Id:		Name:		Type:	Group Name:		Search			
ID	LOGIN_ID	NAME	USER_TYPE	GROUP_ID	GROUP_NAME	CREATED_AT	UPDATED_AT	DELETED_AT	DELETED	USER MODIFY
<input type="checkbox"/>	3	test	RESEARCHER	1	administrators	2017-04-28T01:53:28.396955Z	-	-	N	

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
User Modify


Login ID	test
Name	test
Password	Password
User Type	RESEARCHER
User Group	administrators

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
ID	LOGIN_ID	NAME	USER_TYPE	GROUP_ID	GROUP_NAME	CREATED_AT	UPDATED_AT	DELETED_AT	DELETED	USER MODIFY
<input checked="" type="checkbox"/>	3	test	RESEARCHER	1	administrators	2017-04-28T01:53:28.396955Z	-	-	N	
<input type="checkbox"/>	2	controller	ADMIN	1	administrators	2016-05-19T00:55:30.748446Z	-	-	N	
<input type="checkbox"/>	1	admin	ADMIN	1	administrators	2016-05-19T00:55:30.748446Z	-	-	N	

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
ID	LOGIN_ID	NAME	USER_TYPE	GROUP_ID	GROUP_NAME	CREATED_AT	UPDATED_AT	DELETED_AT	DELETED	USER MODIFY
<input type="checkbox"/>	3	test	RESEARCHER	1	administrators	2017-04-28T01:53:28.396955Z	-	-	N	
<input type="checkbox"/>	2	controller	ADMIN	1	administrators	2016-05-19T00:55:30.748446Z	-	-	N	
<input type="checkbox"/>	1	admin	ADMIN	1	administrators	2016-05-19T00:55:30.748446Z	-	-	N	

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The screenshot shows the "Server Monitoring" page of the NGeneAnalySys Manager. The left sidebar has a red box around the "Server Monitoring" option, which is selected. A red arrow points from this sidebar to the main content area. The main content area contains four sections: 1. CPU usage graph (CPU: Intel(R) Xeon(R) CPU E5-2690 0 @ 2.90GHz core: 1(12 thread)) with a red box around it. 2. Memory usage graph (Total Memory: 15.61GB) with a red box around it. 3. Network traffic section showing Send (9.302 kbyte) and Receive (9.302 kbyte) with a red box around it. 4. Disk usage section showing Disk 1 (total: 18.79GB, 71.4% used) and Disk 2 (total: 3031.63GB, 26.2% used) each with a donut chart showing used (red) and free (green) space.

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The screenshot shows the Analysis Management page of the NGeneAnalySys Manager. At the top left, there is a vertical menu with options: User Management, Server Monitoring, Analysis Management (which is highlighted with a red box), and Log List. Below the menu is a search bar and a table titled "Analysis Job". The table has columns for Run ID, Run Name, Run Status, Samples, User ID, Name, Request Date, Created At, and Job Delete. All rows in the table show a status of "COMPLETE". There are 8 rows in total. A red box highlights the entire table area.

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
This screenshot of the Analysis Management page includes two numbered callouts. Callout 1, located at the top left, points to the search bar with the number "1". Callout 2, located in the middle left, points to the column headers of the table with the number "2". The table structure is identical to the one in the previous screenshot, showing completed analysis jobs.

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
The screenshot shows the 'Analysis Management' section of the NGeneAnalySys Manager. A search bar at the top includes fields for Run Name, Run Status (set to ALL), Request Date, User ID, User Login ID, and User Name, followed by a 'Search' button. Below the search bar is a table titled 'Analysis Job' with the following columns: Run ID, Run Name, Run Status, Samples, User ID, Name, Request Date, Created At, and Job Delete. The table lists eight entries, each with a 'Job Delete' button. The first entry (Run ID 8) has its 'Job Delete' button highlighted with a red box.

Run ID	Run Name	Run Status	Samples	User ID, Name	Request Date	Created At	Job Delete
8	2017-04-28-17-39	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오후 5:39:35	Job Delete
7	2017-04-28-17-24	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오후 5:24:16	Job Delete
6	2017-04-28-17-13	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오후 5:13:36	Job Delete
5	2017-04-28-11-21	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오전 11:21:49	Job Delete
4	2017-04-28-11-02	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오전 11:02:31	Job Delete
3	2017-04-28-11-00	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오전 11:00:35	Job Delete
2	2017-04-28-10-58	FAIL	1	3, test, test	2017-04-28	2017. 4. 28. 오전 10:58:23	Job Delete
1	2017-04-28-10-57	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오전 10:57:11	Job Delete

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This screenshot is identical to the one above, showing the 'Analysis Job' table. The 'Job Delete' button for the first row (Run ID 8) is highlighted with a red box.

Run ID	Run Name	Run Status	Samples	User ID, Name	Request Date	Created At	Job Delete
8	2017-04-28-17-39	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오후 5:39:35	Job Delete
7	2017-04-28-17-24	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오후 5:24:16	Job Delete
6	2017-04-28-17-13	COMPLETE	1	3, test, test	2017-04-28	2017. 4. 28. 오후 5:13:36	Job Delete

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The screenshot shows the 'Log List' page of the NGeneAnalysSys Manager. At the top, there are several search filters: 'Log Type' (set to 'ALL'), 'Log Msg' (empty), 'User ID' (empty), 'User Login ID' (empty), 'User Type' (set to 'ALL'), and 'User Name' (empty). A red box highlights the 'Search' button and these filter fields. Below the filters is a table with 10 rows of log entries. The columns are: Log ID, Log Type, Log Msg, User ID, User Login ID, User Type, Name, and Created At. Each row contains a log entry with 'Log Type' as 'LI' and 'Log Msg' as 'Token issue'. The 'User Type' column consistently shows 'RESEARCHER'. The 'Name' column shows 'test'. The 'Created At' column shows various dates and times from April 28, 2017, such as '2017. 4. 28. 오후 5:59:23' and '2017. 4. 28. 오후 5:44:23'. At the bottom right of the table, there is a navigation bar with numbers 1 through 9 and arrows for pagination.

Log ID	Log Type	Log Msg	User ID	User Login ID	User Type	Name	Created At
88	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:59:23
87	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:54:23
86	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:44:23
85	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:44:23
84	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:39:22
83	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:35:35
82	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:30:35
81	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:29:34
80	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:21:58
79	LI	Token issue	3	test	RESEARCHER	test	2017. 4. 28. 오후 5:16:58

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Code	Country	Code	Country
KOR	South Korea	LBR	Liberia
USA	United States	LBY	Libya
GBR	United Kingdom	LIE	Liechtenstein
CAN	Canada	LTU	Lithuania
JPN	Japan	LUX	Luxembourg
CHN	China	MAC	Macao
AFG	Afghanistan	MKD	Macedonia
ALB	Albania	MDG	Madagascar
DZA	Algeria	MWI	Malawi
ASM	American Samoa	MYS	Malaysia
AND	Andorra	MDV	Maldives
AGO	Angola	MLI	Mali
AIA	Anguilla	MLT	Malta
ATA	Antarctica	MHL	Marshall Islands
ATG	Antigua and Barbuda	MRT	Mauritania
ARG	Argentina	MUS	Mauritius
ARM	Armenia	MYT	Mayotte
ABW	Aruba	MEX	Mexico
AUS	Australia	FSM	Micronesia
AUT	Austria	MDA	Moldova
AZE	Azerbaijan	MCO	Monaco
BHS	Bahamas	MNG	Mongolia
BHR	Bahrain	MNE	Montenegro
BGD	Bangladesh	MSR	Montserrat
BRB	Barbados	MAR	Morocco
BLR	Belarus	MOZ	Mozambique
BEL	Belgium	MMR	Myanmar
BLZ	Belize	NAM	Namibia
BEN	Benin	NRU	Nauru
BMU	Bermuda	NPL	Nepal
BTN	Bhutan	NLD	Netherlands
BOL	Bolivia	ANT	Netherlands Antilles

BIH	Bosnia and Herzegovina	NCL	New Caledonia
BWA	Botswana	NZL	New Zealand
BRA	Brazil	NIC	Nicaragua
IOT	British Indian Ocean Territory	NER	Niger
VGB	British Virgin Islands	NGA	Nigeria
BRN	Brunei	NIU	Niue
BGR	Bulgaria	PRK	North Korea
BFA	Burkina Faso	MNP	Northern Mariana Islands
BDI	Burundi	NOR	Norway
KHM	Cambodia	OMN	Oman
CMR	Cameroon	PAK	Pakistan
CPV	Cape Verde	PLW	Palau
CYM	Cayman Islands	PSE	Palestine
CAF	Central African Republic	PAN	Panama
TCD	Chad	PNG	Papua New Guinea
CHL	Chile	PRY	Paraguay
CXR	Christmas Island	PER	Peru
CCK	Cocos Islands	PHL	Philippines
COL	Colombia	PCN	Pitcairn
COM	Comoros	POL	Poland
COK	Cook Islands	PRT	Portugal
CRI	Costa Rica	PRI	Puerto Rico
HRV	Croatia	QAT	Qatar
CUB	Cuba	COG	Republic of the Congo
CUW	Curacao	REU	Reunion
CYP	Cyprus	ROU	Romania
CZE	Czech Republic	RUS	Russia
COD	Democratic Republic of the Congo	RWA	Rwanda
DNK	Denmark	BLM	Saint Barthelemy
DJI	Djibouti	SHN	Saint Helena
DMA	Dominica	KNA	Saint Kitts and Nevis
DOM	Dominican Republic	LCA	Saint Lucia
TLS	East Timor	MAF	Saint Martin
ECU	Ecuador	SPM	Saint Pierre and Miquelon
EGY	Egypt	VCT	Saint Vincent and the Grenadines

SLV	El Salvador	WSM	Samoa
GNQ	Equatorial Guinea	SMR	San Marino
ERI	Eritrea	STP	Sao Tome and Principe
EST	Estonia	SAU	Saudi Arabia
ETH	Ethiopia	SEN	Senegal
FLK	Falkland Islands	SRB	Serbia
FRO	Faroe Islands	SYC	Seychelles
FJI	Fiji	SLE	Sierra Leone
FIN	Finland	SGP	Singapore
FRA	France	SXM	Sint Maarten
PYF	French Polynesia	SVK	Slovakia
GAB	Gabon	SVN	Slovenia
GMB	Gambia	SLB	Solomon Islands
GEO	Georgia	SOM	Somalia
DEU	Germany	ZAF	South Africa
GHA	Ghana	SSD	South Sudan
GIB	Gibraltar	ESP	Spain
GRC	Greece	LKA	Sri Lanka
GRL	Greenland	SDN	Sudan
GRD	Grenada	SUR	Suriname
GUM	Guam	SJM	Svalbard and Jan Mayen
GTM	Guatemala	SWZ	Swaziland
GGY	Guernsey	SWE	Sweden
GIN	Guinea	CHE	Switzerland
GNB	Guinea-Bissau	SYR	Syria
GUY	Guyana	TWN	Taiwan
HTI	Haiti	TJK	Tajikistan
HND	Honduras	TZA	Tanzania
HKG	Hong Kong	THA	Thailand
HUN	Hungary	TGO	Togo
ISL	Iceland	TKL	Tokelau
IND	India	TON	Tonga
IDN	Indonesia	TTO	Trinidad and Tobago
IRN	Iran	TUN	Tunisia
IRQ	Iraq	TUR	Turkey
IRL	Ireland	TKM	Turkmenistan
IMN	Isle of Man	TCA	Turks and Caicos Islands

ISR	Israel	TUV	Tuvalu
ITA	Italy	VIR	U.S. Virgin Islands
CIV	Ivory Coast	UGA	Uganda
JAM	Jamaica	UKR	Ukraine
JEY	Jersey	ARE	United Arab Emirates
JOR	Jordan	URY	Uruguay
KAZ	Kazakhstan	UZB	Uzbekistan
KEN	Kenya	VUT	Vanuatu
KIR	Kiribati	VAT	Vatican
XKX	Kosovo	VEN	Venezuela
KWT	Kuwait	VNM	Vietnam
KGZ	Kyrgyzstan	WLF	Wallis and Futuna
LAO	Laos	ESH	Western Sahara
LVA	Latvia	YEM	Yemen
LBN	Lebanon	ZMB	Zambia
LSO	Lesotho	ZWE	Zimbabwe

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Code	Ethnic Origin
ARAB_LEAGUE	Arab League
EUROPE	Europe
AUSTRALIAN_INDIGENOUS	Australian (indigenous)
INDIGENOUS_OCEANIAN	Indigenous Oceanian
EUROPEANS_IN_OCEANIA	Europeans in Oceania
INDIGENOUS	Indigenous
CANADA	Canada
MEXICO	Mexico
UNITED_STATES	United States
CENTRAL_AMERICA	Central America
SOUTH_AMERICA	South America
ASIA_CENTRAL	Asia-Central
ASIA_EAST	Asia-East
ASIA_NORTHERN	Asia-Northern
ASIA_SOUTH	Asia-South
ASIA_SOUTHEAST	Asia-Southeast
ASIA_WEST	Asia-West

PATHOGENICITY CLASS

Pathogenic

- Nonsense mutations (stop_gained) with 1000 genomes and ESP6500 <= 0.01.
- Frameshift (frameshift_variant) with 1000 genomes and ESP6500 <= 0.01.
- Splice acceptor (splice_acceptor_variant) sites (-2 and -1) and/or splice donor (splice_donor_variant) sites (+1 and +2).
- Loss of natural initiation of transcription (start codon, start_lost) or termination (stop codon, stop_lost).
- In frame INDELs (inframe_insertion, inframe_deletion) of more than 3 nucleotides.
- Mutation reported as pathogenic in BIC or ClinVar.

Likely Pathogenic

- Missense (missense_variant) mutation with 1000 genomes and ESP6500 < 0.01.
- In frame INDELs of 3 nucleotides.
- Missense or synonymous with at least 2 values of predictive score (SIFT, PolyPhen2, MutationTaster) > 0.9.

Likely Benign

- None.

Benign

- Missense, synonymous or intronic variant with 1000 genomes or ESP6500 > 0.05 with 2 predictive.

Uncertain Significance

- Mutation no classified in pathogenic, likely pathogenic and benign.

PATHOGENICITY CLASS - ClinVar

Clinical significance on ClinVar submissions (SCV)

Clinical Significance Value	Guidance for Use in ClinVar SCV Records
Benign	As recommended by ACMG/AMP for variants interpreted for Mendelian disorders.
Likely benign	As recommended by ACMG/AMP for variants interpreted for Mendelian disorders.
Uncertain significance	As recommended by ACMG/AMP for variants interpreted for Mendelian disorders.
Likely pathogenic	As recommended by ACMG/AMP for variants interpreted for Mendelian disorders.
Pathogenic	As recommended by ACMG/AMP for variants interpreted for Mendelian disorders.
drug response	A general term for a variant that affects a drug response, not a disease. We anticipate adding more specific drug response terms based on a recommendation by CPIC .
association	For variants identified in a GWAS study and further interpreted for their clinical significance.
risk factor	For variants that are interpreted not to cause a disorder but to increase the risk.
protective	For variants that decrease the risk of a disorder, including infections.
Affects	For variants that cause a non-disease phenotype, such as lactose intolerance.
conflicting data from submitters	Only for submissions from a consortium, where groups within the consortium have conflicting interpretations of a variant but provide a single submission to ClinVar.
other	If ClinVar does not have the appropriate term for your submission, we ask that you submit "other" as clinical significance and contact us to discuss if there are other terms we should add.
not provided	For submissions without an interpretation of clinical significance. The primary goal of ClinVar is to archive reports of clinical significance of variants. Therefore, submissions with a clinical significance of "not provided" should be limited to: <ul style="list-style-type: none"> • "literature only" submissions that report a publication about the variant, without interpreting the clinical significance.

- “research” submissions that provide functional significance (e.g. undetectable protein level) but no interpretation of clinical significance.
- “clinical testing” or “phenotyping only” submissions from clinics or physicians that provide additional information about individuals with the variant, such as observed phenotypes, but do not interpret the clinical significance.

(<https://www.ncbi.nlm.nih.gov/clinvar/docs/clinsig/>)

Clinical significance on ClinVar aggregate records (RCV)

Combination of Values from Different Submitters	Reported as
Pathogenic and Likely pathogenic	Pathogenic/Likely Pathogenic
(Pathogenic or Likely pathogenic or Benign or Likely benign) AND Uncertain significance	Conflicting interpretations of pathogenicity
(Pathogenic or Likely pathogenic) AND (Benign or Likely benign)	Conflicting interpretations of pathogenicity
Benign and Likely benign	Benign/Likely benign
Any single ACMG value AND any non-ACMG value, e.g. Uncertain significance and risk factor	Uncertain significance, risk factor
Conflicting ACMG values AND any non-ACMG value, e.g. Pathogenic and Uncertain significance and risk factor	Conflicting interpretations of pathogenicity, risk factor
No ACMG value AND multiple non-ACMG values, e.g. Drug response and risk factor	drug response, risk factor
A consortium does not agree on which ACMG assessment is valid	conflicting data from submitters

(<https://www.ncbi.nlm.nih.gov/clinvar/docs/clinsig/>)

PATHOGENICITY CLASS - BIC

Class	Clinical Significance Value	Description
Class 1	Not pathogenic low clinical significance	There is significant evidence against this variant being a dominant high-risk pathogenic mutation.
Class 2	Likely Notice pathogenic little clinical significance	There is evidence against this variant being a dominant high-risk pathogenic mutation.
Class 3	Uncertain	There is insufficient evidence to place this variant in Class 1,2,4 or 5.
Class 4	Likely Pathogenic	There is evidence that this variant is a dominant high-risk pathogenic mutation.
Class 5	Pathogenic	There is significant evidence to suggest that this variant is a dominant high-risk pathogenic mutation.

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번호	필드명	Defalut	설명
사용자 설정 항목			
1	Prediction	Y	Pathogenic class using ACMG ISV guideline.
2	Pathogenicity	Y	User defined Pathogenic class.
3	Filter	Y	Pre-filtered variant in BRCA1/2.
4	Warn	Y	Low confidence variants.
5	Report	Y	Reported variant.
6	False	Y	Set variant to false.
7	ID	N	Variant number.
변이 및 유전자 annotation 정보			
8	Type	Y	Variant type SNP: Single-Nucleotide Polymorphism (R='A', S='C'). INS: Insertion (R='A', S='AGT'). MNP: Multiple-nucleotide polymorphism (R='ATA', S='GTC'). MIXED: Multiple-nucleotide and an Indel (R='ATA', S='GTCAGT').
9	Cod.Cons	Y	
10	Gene	Y	Gene Symbol (HUGO Gene Nomenclature).
11	Strand	Y	Whether the gene is located on the plus or minus strand.
12	Transcript	Y	Transcript ID.
13	NT change	Y	The variant coordinates relative to the coding DNA (HGVS format: A of start codon ATG is #1).
14	AA change	Y	Change amino acid.
15	NT change(BIC)	Y	The variant coordinates relative to the coding DNA (BRCA format-BRCA1: A of start codon ATG is #120, BRCA2: A of start codon ATG is #229).
16	ref.genome	N	Reference genome build.
변이 위치 정보			
17	Chr	Y	The chromosome numbers.
18	Ref	Y	Reference allele.
19	Alt	Y	Alternative allele.
20	Zygosity	Y	The words homozygous, heterozygous, and hemizygous are used to describe the genotype of a diploid organism at a single locus on the DNA.
21	g.pos	Y	The variant coordinate on the reference genome.
22	g.pos.end	N	The variant end coordinate on the reference genome.
23	Exon	Y	Number of exon in which the variant is located exon/total

			exon.
24	Exon(BIC)	Y	Exon number; Reference sequences: BRCA1 GenBank U14680; BRCA2 GenBank U43746.
Allele 정보			
25	ref.num	Y	Number of reads supporting the reference allele.
26	alt.num	Y	Number of reads supporting the alternative allele.
27	depth	Y	Total coverage of the variant.
28	fraction	Y	The proportion of the reads supporting the variant count.
Population 빈도			
29	dbSNP	N	dbSNP RS ID.
30	1KG	Y	1000 Genomes Project population frequency.
31	ExAC	Y	Exome Aggregation Consortium(ExAC) population frequency.
32	Esp6500	Y	NHLBI Exome Sequencing Project(ESP) population frequency.
33	Korean	Y	Korean Exome Project frequency.
데이터베이스 정보			
34	ClinVar.Acc	Y	ClinVar RCV accessions.
35	ClinVar.disease	N	ClinVar disease name.
36	ClinVar.Class	Y	Clinical classification.
37	BIC.Cat	N	Mutation type in BIC. 3UTR: 3'UTR. 5TR: 5' UTR. F: Frameshift. IFI: In frame insertion. IFD: In frame deletion. IVS: Intervening sequence. M: Missense. N: Nonsense. P: Polymorphism. S: Splice. Syn: Synonymous. UV: Unclassified variant.
38	BIC.Importance	N	Clinical importance in BIC.
39	BIC.Class	Y	Clinical classification using 5 classes in BIC. 0: pending. 1: Not pathogenic/low clinical significance. 2: Likely not pathogenic/little clinical significance. 3: Uncertain. 4: Likely pathogenic. 5: Pathogenic.
40	Experiment	N	Experiment type {somatic, germline}.

41	BIC.Designation	Y	Designation of described mutation according to BIC nomenclature guidelines.
42	BIC.NT	Y	Nucleotide # in the transcript (cDNA) at which mutation occurs; Reference sequences: BRCA1 GenBank U14680; BRCA2 GenBank U43746.
43	KOHBRA.patient	Y	High risk breast cancer patient in Korea (KOHBRA).
44	KOHBRA.frequency	Y	High risk breast cancer frequency in Korea (KOHBRA).
45	Be.clinvar.update	N	Date the variant was last updated in ClinVar.
46	Be.clinvar.origin	N	Variants are classified as either germline or somatic, depending on how they are acquired. Germline variants are genetic changes that we inherit from our parents. Somatic variants are DNA changes that we acquire over our lifetime, often through exposure to pollutants, toxins, radiation and other carcinogens.
47	Be.clinvar.method	N	Method used to collect the data that supports the assertion of clinical significance. Allowed values: case-control, clinical testing, literature only, reference population, research.
48	Be.BIC.category	N	Mutation type in BIC.
49	Be.BIC.nationality	N	Nationality of the patient associated with the variant. From BIC.
50	Be.BIC.ethnic	N	Ethnic of BIC.
51	Be.BIC.patho	N	BIC pathogenicity.
52	Be.transcript	N	Reference transcript.
53	Be.nt	N	Nucleotide change.
54	Be.gene	N	Reference gene.
55	Be.enigma.condition	N	Human-readable condition, describing the biological impact of the variation.
56	enigma	N	
57	Be.enigma.update	N	The date on which the clinical significance of the variant was last evaluated by ENIGMA.
58	Be.exLOVD.class	N	IRAC Classification from exLOVD.
59	Be.clinvar.patho	N	ClinVar pathogenicity.
60	Be.enigma.patho	Y	The Clinical Significance columns indicates whether expert curators have determined if the variant is pathogenic or benign. What do these classifications mean? Pathogenic: variants confer an increased risk of disease. Likely pathogenic: variants have good evidence to support an association with disease risk. Likely benign: variants have good evidence to support no association with disease risk. Benign: variants are not associated with any markedly increased risk of disease. Variants of uncertain significance (VUS): are those for which the evidence of disease risk is

			not clear yet, sometimes because there is not yet enough evidence to classify them as either pathogenic or benign.
61	polyphen2	Y	PolyPhen-2 score ranges from 0.0 (tolerated) to 1.0 (deleterious).
62	sift	Y	SIFT score ranges from 0.0 (deleterious) to 1.0 (tolerated).
63	mutationtaster	Y	MutationTaster score.
64	Left.seq	N	Left 22bp sequence in variant.
65	Right.seq	N	Right 22bp sequence in variant.

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- 1) Lynch HT, Snyder CL, Lynch JF, Riley BD, Rubinstein WS. Hereditary breast-ovarian cancer at the bedside: role of the medical oncologist. *Journal of clinical oncology*. 2003;21(4):740-53.
- 2) Robson M, Offit K. Clinical practice. Management of an inherited predisposition to breast cancer. *The New England journal of medicine*. 2007;357(2):154-62.
- 3) Pal T, Permuth-Wey J, Betts JA, et al. BRCA1 and BRCA2 mutations account for a large proportion of ovarian carcinoma cases. *Cancer*. 2005;104(12):2807-16.
- 4) Kang E, Seong MW, Park SK, et al. The prevalence and spectrum of BRCA1 and BRCA2 mutations in Korean population: recent update of the Korean Hereditary Breast Cancer (KOHBRA) study. *Breast Cancer Res Treat*. 2015;151:157-68.
- 5) Moyer VA. Risk assessment, genetic counseling, and genetic testing for BRCA-related cancer in women: U.S. Preventive Services Task Force recommendation statement. *Annals of internal medicine*. 2014;160(4):271-81.
- 6) Gadzicki D, Evans DG, Harris H, et al. Genetic testing for familial/hereditary breast cancer-comparison of guidelines and recommendations from the UK, France, the Netherlands and Germany. *Journal of community genetics*. 2011;2(2):53-69.
- 7) American Society of Clinical Oncology policy statement update: genetic testing for cancer susceptibility. *Journal of clinical oncology*. 2003;21(12):2397-406.
- 8) Rechards S, Aziz N, Bale S, et al. Standards and guidelines for the interpretation of sequence variants: a joint consensus recommendation of the American College of Medical Genetics and Genomics and the Association for Molecular Pathology. *Genetics in Medicine*. 2015;17(5):405-24.
- 9) CPIC Term Standardization for Clinical Pharmacogenetic Test Results Project.
<https://www.pharmgkb.org/page/cpicTermProject>
- 10) Representation of clinical significance in ClinVar and other variation resources at NCBI.
<https://www.ncbi.nlm.nih.gov/clinvar/docs/clinsig/>

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Copyright © 2016 by NGeneBio Co., Ltd. All Rights Reserved.